

Product datasheet for MC209321

Rp1 (NM_001195662) Mouse Untagged Clone

Product data:

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product Type:	Expression Plasmids
Product Name:	Rp1 (NM_001195662) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Rp1
Synonyms:	Dcdc3; Gm38717; mG145; O; Orp1; Rp; Rp1h
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>MC209321 representing NM_001195662 Red=Cloning site Blue=ORF Orange=Stop codon
	TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC GCC <mark>GCGATCGC</mark> C
	ATGCTTGAGAAATTGCAGGTCTCACCCAAAATGAGTGACACACCTTCTACTAGTTTCTCCATGATTCATC TGACTTCTGAAGGTCAAGTTCCTTCCCCTCGCCATTCAAATATCACTCATCCTGTAGTGGCTAAACGCAT CAGTTTCTATAAGAGTGGAGACCCACAGTTTGGCGGCGTTCGGGTGGTGACCCCTCGTTCCTTTAAG ACTTTTGACGCTCTGCTGGACAGTTTATCCAGGAAGGTACCCCTGCCCTTTGGGGTAAGGAACATCAGCA CGCCCCGTGGACGACACAGCATCACCAGGCTGGAGGAGGTACCCCTGCCCTTTGGGGTAAGGAACATCAGCA CGCCCCGTGGACGACACAGCATCACCAGGCTGGAGGAGGTAGGGACGGCAAGTCTTATGTGTGCTCCCA CAATAAGAAGGTGCTGCCAGTTGACCTGGACAAGGCCCGCAGGCGCCCTCGGCCCTGGCCTGAGTAGTCGC TCCATAAGCACGCATGTGCAGCTCTGTCCTGCAACTGCCAATATGTCCACCATGGCACCTGGCATGCTCC GTGCCCCAAGGAGGCTCGTGGTCTTCCGGAATGGTGACCCGAAGAATAAGCATGTGGTCCTTCTTAGCCG GAGGATCACTCAGAGCTTCGAAGCCTTTCTGCAGTACCTGACACAGGTTATGCAGTGTGCTCCTGTGGCCAAG CTGTATGCCACAGATGGAAGAAAAGTTCCTAGTCTCCAGGCAGTGATACTCAGCTCTGGAGCTGTGGTGG CAGCTGGAAGGGAGCCATTTAAACCAGGAAATTACGACATACAAAAGTACTTGCTTCCTGCCAAGTTACC AGGAATCTCTCATCGTGTGCACCAGAAGGGAAAAGCCAAGATAGAAAAAGAAAG
Restriction Sites:	Sgfl-Mlul
ACCN:	NM_001195662
Insert Size:	858 bp



This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2024 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US

	Rp1 (NM_001195662) Mouse Untagged Clone – MC209321
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution M	 ethod: 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM 001195662.1, NP 001182591.1</u>
RefSeq Size:	3047 bp
RefSeq ORF:	858 bp
Locus ID:	19888
Cytogenetics:	1 1.65 cM
Gene Summary:	This gene encodes a member of the doublecortin family. The protein encoded by this gene contains two doublecortin domains, which bind microtubules and regulate microtubule polymerization. The encoded protein is a photoreceptor microtubule-associated protein and is required for correct stacking of outer segment disc. This protein and the RP1L1 protein, another retinal-specific protein, play essential and synergistic roles in affecting photosensitivity and outer segment morphogenesis of rod photoreceptors. Because of its response to in vivo retinal oxygen levels, this protein was initially named ORP1 (oxygen- regulated protein-1). This protein was subsequently designated RP1 (retinitis pigmentosa 1) when it was found that mutations in this gene cause autosomal dominant retinitis pigmentosa. Mutations in this gene also cause autosomal recessive retinitis pigmentosa. [provided by RefSeq, Jun 2019] Transcript Variant: This variant (2) has an alternate 5' exon, which results in an upstream in- frame AUG start codon, and an alternate 3' exon, as compared to variant 1. The resulting isoform (2) has a longer N-terminus, but has a much shorter and distinct C-terminus, as compared to isoform 1. The isoform 2 retains two doublecortin domains.

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2024 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US