

Product datasheet for **MC209361**

Ccl1 (NM_011329) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Ccl1 (NM_011329) Mouse Untagged Clone
Tag: Tag Free
Symbol: Ccl1
Synonyms: BF534335; I-309; P500; Scya1; Tca-3; Tca3
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Fully Sequenced ORF: >MC209361 representing NM_011329
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGAAACCCACTGCCATGGCACTGATGTGCCTGCTGCTGGCTGCCGTGTGGATACAGGATGTTGACAGCA
AGAGCATGCTTACGGTCTCCAATAGCTGCTGCTGAACACCTTGAAGAAAGAGCTCCCTGAAGTTTAT
CCAGTGTACAGAAAGATGGGCTCCTCCTGTCTGATCCCCAGCTGTGGTATTCAGGCTGAACAAAGGT
AGAGAAAGCTGCGCCTCAACTAACAAAACGTGGGTTCAAATCACCTGAAGAAGGTGAACCCCTGCT**AA**

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Chromatograms: https://cdn.origene.com/chromatograms/ja3096_a03.zip

Restriction Sites: SgfI-MluI

ACCN: NM_011329

Insert Size: 279 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_011329.3</u> , <u>NP_035459.1</u>
RefSeq Size:	505 bp
RefSeq ORF:	279 bp
Locus ID:	20290
UniProt ID:	<u>P10146</u>
Cytogenetics:	11 49.99 cM
Gene Summary:	Cytokine that is chemotactic for neutrophils.[UniProtKB/Swiss-Prot Function]