

## Product datasheet for **MC209385**

### Sftpc (NM\_011359) Mouse Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Sftpc (NM_011359) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Sftpc
Synonyms:	Bricd6; pro-SpC; Sftp-2; Sftp2; SP-C; SP5; SPC
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>MC209385 representing NM_011359 Red=Cloning site Blue=ORF

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGC**C

ATGGACATGAGTAGCAAAGAGGTCCTGATGGAGAGTCCACCGGATTACTCGGCAGGTCCCAGGAGCCAGT  
TCCGCATCCCCTGCTGTCCCGTGCACCTCAAACGCCTTCTCATCGTGGTTGTGGTGGTGGTCTCGTTGT  
CGTGGTGATTGTAGGGGCTCTGCTCATGGGCTCCACATGAGTCAAAAACATACTGAGATGGTCCTTGAG  
ATGAGCATCGGAGCACCGAAACTCAGAAACGCCTAGCCCGAGTGAGCGAGCAGACACCATCGCTACCT  
TTCCATCGGCTCCACTGGCATCGTTGTGTATGACTACCAGCGGCTCCTGACGGCTATAAGCCAGCTCC  
AGGAACCTACTGCTACATCATGAAGATGGCTCCAGAGAGCATCCCTAGTCTTGAGCTTTTCGCTAGAAAA  
CTCCAGAACTTCCAGGCCAAGCCCTCCACACCCACCTTAAGCTGGGCCAGGAGGAAGGGCATGATACTG  
GTTCCGAGTCCGATTCTTCCGGGAGAGACCTGGCTTTCCTAGGCCTTGCTGTGAGCACCTGTGTGGAGA  
GCTACCACTCTACTATATCTAG

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

Chromatograms:	<a href="https://cdn.origene.com/chromatograms/ja1603_h10.zip">https://cdn.origene.com/chromatograms/ja1603_h10.zip</a>
Restriction Sites:	Sgfl-Mlul
ACCN:	NM_011359
Insert Size:	582 bp



[View online »](#)

<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">BC061137</a></u> , <u><a href="#">AAH61137</a></u>
<b>RefSeq Size:</b>	819 bp
<b>RefSeq ORF:</b>	582 bp
<b>Locus ID:</b>	20389
<b>Cytogenetics:</b>	14 36.32 cM