

Product datasheet for **MC209404**

Cox7a2l (NM_001159529) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Cox7a2l (NM_001159529) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Cox7a2l
Synonyms:	COX7AR; COX7RP; EB1; SIG-81; SIG81; Silg81
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin
Fully Sequenced ORF:	>MC209404 representing NM_001159529 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGTACTACAAGTTTAGCAGTTTCACGCAGAAGTTGGCTGGAGCTTGGGCTTCGGAAGCCTACACCCCGC
AGGGGTTAAAGCCAGTTCCACAGAAGCACCATCATATTTGCCACCAACCAAACACTGACCTCCAG
TGTGACAGCATATGATTATTCTGGGAAGAACAAGTTCCAGAGCTGCAGAAGTTTTCCAGCCTCCTTTT
CTTAAGGAGAAGCTTTACTTCAGTTTGGTCCGTTGGCTAATAATGGAACCTTTGTACAGAAGGCTGATG
GTTTCCACCTGAAACGAGGCCTTCAGACCAAATGCTTTACCGGACCACCATGGCTCTGACTGGGAGG
GACCATCTACTGCCTGATCGCCCTCTACATGGCCTCGCAGCCAGAAACAAATGA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	Sgfl-MluI
ACCN:	NM_001159529
Insert Size:	405 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).



[View online »](#)

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001159529.1](#), [NP_001153001.1](#)

RefSeq Size: 1169 bp

RefSeq ORF: 405 bp

Locus ID: 20463

Cytogenetics: 17 E4

Gene Summary: Involved in the regulation of oxidative phosphorylation and energy metabolism (PubMed:23857330, PubMed:23812712). Necessary for the assembly of mitochondrial respiratory supercomplex (PubMed:23857330, PubMed:23812712).[UniProtKB/Swiss-Prot Function]
Transcript Variant: This variant (1) represents the longer transcript and encodes the longer isoform (1).