

Product datasheet for MC209423

Sstr4 (NM_009219) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Sstr4 (NM_009219) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Sstr4
Synonyms:	Smstr4; sst4
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>MC209423 representing NM_009219 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCTGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGC**

ATGAACGCGCCAGCAACTCTGCCTCCGGGGTTCGAGGACACCACTGGACCCCTGGGATCAACGCCAGCT
GGGCTCCGGACGAGGAGGAGGATGCGATGGGGTCCGACGGCACAGGGACAGCGGGCATGGTCACTAT
CCAGTGCACTATGCGCTCGTGTCTGGTGGGCTGGTGGGAAACGCCCTGGTCATCTTCGTGATCCTA
CGCTATGCCAAGATGAAGACAGCCACCAACATCTACCTACTCAACCTGGCCGTCGCCGATGAGCTCTTCA
TGCTCAGCGTGCCATTCTGTGGCCTCGGCGGCTGCCCTGCGCCACTGGCCGTTTCGGGCGGTGCTGTGTCTG
CGCAGTGCTTAGCGTGGACGGCCTGAACATGTTCACTAGTGTCTTCTGCCTCACCGTGCTCAGCGTGGAC
CGCTATGTGGCTGTTGTGCACCCCTCTGCGCGCCGCCACCTACCGCGGCCAGCGTGGCCAAGCTAATCA
ACCTGGGAGTGTGGCTAGCATCCTTGCTGGTCAACCTGCCATCGCAGTCTTCGCTGACACCAAGGCCAGC
TCGTGGGGGCGAGGCCGTGGCTTGCAACCTGCACTGGCCTCACCGGCCCTGGTCTGCGGTCTTTGTGATC
TATACTTTTTTGTGGGCTTCTACTCCCGTTCTGGCCATCGGATTATGCTACCTGCTTATTGTGGCA
AGATGCGCGCTGTGGCCTGCGGGCTGGCTGGCAACAACGGAGGCGCTCAGAGAAGAAGTCACTAGGCT
CGTGCTAATGGTGGTGACCGTCTTTGTGCTATGCTGGATGCCTTTCTATGTGGTGCAGCTTCTGAACCTG
TTTGTCACCAAGCTCGATGCCACTGTCAACCATGTGTCCCTCATCCTGAGCTATGCCAAGCTGTGCCA
ATCCCATTCTCTATGGCTTCTCTCTGACAACCTCCGGCGCTTTCCAGCGGGTCTGTGCCTGCGCTG
CTGTCTCTGAAACAACCTGGAGGTGCTGAGGAAGAGCCCTGGACTACTATGCCACTGTCTCAAAGC
AGAGGGGGTGAGGATGCATATGCCCTCTCTGCCCTGCCAGCAGGAGCCGTGCAAGCAGAACCTGGCT
GCAAGCAAGTCCCTTTCACCAAGACCACTACTTCTGA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Chromatograms:	https://cdn.origene.com/chromatograms/ja3041_f01.zip
Restriction Sites:	Sgfl-MluI
ACCN:	NM_009219
Insert Size:	1158 bp
OTI Disclaimer:	Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery. The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	NM_009219.3 , NP_033245.2
RefSeq Size:	1388 bp
RefSeq ORF:	1158 bp
Locus ID:	20608
UniProt ID:	P49660
Cytogenetics:	2 73.44 cM
Gene Summary:	Receptor for somatostatin-14. The activity of this receptor is mediated by G proteins which inhibits adenylyl cyclase. It is functionally coupled not only to inhibition of adenylate cyclase, but also to activation of both arachidonate release and mitogen-activated protein (MAP) kinase cascade.[UniProtKB/Swiss-Prot Function]