

Product datasheet for **MC209431**

Sox14 (NM_011440) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Sox14 (NM_011440) Mouse Untagged Clone
Tag: Tag Free
Symbol: Sox14
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Fully Sequenced ORF: >MC209431 representing NM_011440
Red=Cloning site **Blue**=ORF **Orange**=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGTCCAAACCTTCAGACCACATCAAGCGGCCCATGAATGCCTTCATGGTATGGTCCCGGGGCCAGCGGC
GCAAGATGGCCAGGAAAATCCCAAGATGCACAACCTCGGAGATCAGCAAACGCCCTGGGCGCCGAATGGAA
GCTCCTGTCCGAGGCAGAGAAGCGGCCCTTACATCGATGAAGCCAAAAGGCTGCGCGCCAACACATGAAG
GAGCACCTGACTACAAGTACCGGCCGCGGCCAAGCCGAAGAACCTGCTCAAGAAGGACAGGTATGTCT
TCCCCCTGCCCTACCTGGGCGACACGGACCCGCTCAAGGCGGCCGGCCTGCCCGTGGGGGCTCCGACGG
CCTCCTGAGCGGCCGAGAAAAGCCCGGCCCTTCTTGCCGCTGCTCGGCACCTTACTCCTGCTGGAC
CCTGCGCAATTTAGTTCCAGTGCCATCCAGAAGATGGGTGAAGTGCCCCACACGTTGGCCACCAGCGCGC
TGCCCTACGCGTCCACCCTGGGCTATCAGAACGGCGCCTTCGGCAGTCTCAGCTGCCCCAGCCAGCACAC
GCACACGCACCCGTCCCCACCAACCCGGGCTACGTGGTGCCCTGTAAGTGTACCGCCTGGTCTGCCTCC
ACCCTGCAGCCCCCTGTCGCTACATCCTTCCAGGCATGACCAAGACTGGCATAGACCCTTATTTCGT
CAGCCCACGCTACGGCCAT**GTA**A

AG**CGGACCG**ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCC
TGGATTACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-RsrII
ACCN: NM_011440
Insert Size: 723 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_011440.1 , NP_035570.1
RefSeq Size:	2075 bp
RefSeq ORF:	723 bp
Locus ID:	20669
UniProt ID:	Q04892
Cytogenetics:	9 52.11 cM
Gene Summary:	Acts as a negative regulator of transcription.[UniProtKB/Swiss-Prot Function]