

## Product datasheet for MC209441

### Serpinb9f (NM\_183197) Mouse Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Serpinb9f (NM\_183197) Mouse Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Serpinb9f  
**Synonyms:** BC052216; NK21; NK21L1; ovalbumin; Spi13  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**Fully Sequenced ORF:** >MC209441 representing NM\_183197  
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGAATACTCTGTCTCAAGCCAATGGCACCTTTGCCATCCATCTTTGAAGGTGCTCTGTCAAGACAACC  
CTTCAAAAAATGTGTATTCTCCTATGAGCATCTTCTGCTCTAGCCATGGTCTCTTGGGTGCAA  
GGGAGACTGTCAGTCCAGATATGTCAGGCACTTCATTTAAACCCAGATGAAGATGTCCATCAGGGCTTC  
CAGTTGCTTCTCCATAACTTGAACAAGCAAAACAACAAAAGTACTGCCTTAGAATGGCCAACAGGCTCT  
TTGTAGAAAATACTTGTGAATTACTTCCAACATTTAAGGAGTCCGTCTTAAATTCTACCACTCTGAGAT  
GGAGCAGCTCTCTTTGCCAAAGCTGCCGAGGAGTCCAGGCAACACATAAACATGTGGGTCTCCAAACAG  
ACTAATGGTAAAATCCAGACTTATTGTCAAAGACTCAGTTAATTCACAGACCAGGCTGATTCTCGCCA  
ATGCCTTATATTTCCATGGAACATGGTGTAAACGTTTTGAAAAAATAGGACCAAGGAAATGCCCTTTAA  
GATAACAAGAAGGAAACAAGACCAGTGCAGATGATGTGGCGGGAAGACACATTATCCATGCCTATGTG  
AAGGAAATCCAGGCACAAGTGTAGTATGCCCTATGAGGCATAGATTTGAATTTTGGTCTACTCC  
CAGATGAGGGTGTGGACATCAGCAAGGTGGAAAACAATCTCACTTTTGAGAAGTTAACAGCCTGGACCAA  
ACCAGAGTTTATGAATCGCACTGAGTTTCATGTTTATCTTCAAAAATTTCAACTACAAGAGGATTATGAT  
ATGAATCCCTGCTACAGCATTGGGAATTTGAAATGTCTTCGATGGAAGCAAGGCTGACTTATCAGGAA  
TGTCTACTAAGGAAAACCTGTGTCTGTCCGAGTTTGTTCACAAGTGTGTAGTAGAGGTCATGAGGAAGG  
AACAGAGGCTGTGCAGCCTCTGCTGTAGAATTTATTTTTATGTTTAGGTCCTGATCCTGAAACTTTCT  
TGTGCTGACCACCCCTTCTTTTCTTCATCATGCACAGCACAACTAACAGCATCCTGTTCTGTGGCAGAT  
TCTCATCTCC**TAA**

**ACGGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** SgfI-MluI



[View online >](#)

<b>ACCN:</b>	NM_183197
<b>Insert Size:</b>	1134 bp
<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">NM_183197.1</a>, <a href="#">NP_899020.1</a></u>
<b>RefSeq Size:</b>	1920 bp
<b>RefSeq ORF:</b>	1134 bp
<b>Locus ID:</b>	20709
<b>Cytogenetics:</b>	13 13.86 cM