

Product datasheet for **MC209462**

Sprr3 (NM_011478) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Sprr3 (NM_011478) Mouse Untagged Clone
Tag: Tag Free
Symbol: Sprr3
Synonyms: SPR3
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Fully Sequenced ORF: >MC209462 representing NM_011478
Red=Cloning site **Blue**=ORF **Orange**=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGAGTTCCTTACCAGCAAAGCAGCCCTTTGTCCCACCTCCTCAGCCTGAACAGCATCAAGTGAAGCAAC
CTTGTCAGCCTCCACCACAAGGAAAATTTGTTCCCATAGCCACCTCTGAGCCATGCCACACAGATGTTCC
ACAACCAGGAAACACCAAGATTCAGAACCATGCAGCACCAAAGTGCCTGAGCCAGGGAACACTGTTGTT
CTTGAGCCAGATTATACTACGATGCCTGGGCCATGTTCTACAAACATCACTGAGCCAGATTATACCACAA
TTCTGGGCCATGCTCTACCAACATCACAGAGCCAGATTATACTACAATTCCTGGGCCATGCTCTACCAA
CATCCCTGGGCCAGATCGTACCGTGGTTCTGGATCATGCTCTACCAACATCACAGAGCCAGATTATACC
ACAATTCCTGGGCCATCCTCTACCAAAATCCCTGACCCAGGCTGTGCCATGGTTCCTGGGCCAGTCCCT
CCAGTACTTCTGAGCCAAGCTCAGAACCATGCTCTATCAATGTTTCGTGAACCTGGCTACATGAATGCCTC
AGAGCCCACCCATGCCAAGGTTCTGATCAGGGCTACACCAAGATACCTGACCAGGGATCCTCTAAGGTA
CCAGAGCCATGCCAGTCCAGGGTTCTGAGGTGTGCCCCCAACTGTCAACCCTGTATCAGCTAAGCAGA
AGACCAAGCAGAAG**TAA**

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_011478
Insert Size: 717 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_011478.2 , NP_035608.1
RefSeq Size:	1222 bp
RefSeq ORF:	717 bp
Locus ID:	20766
UniProt ID:	O09116
Cytogenetics:	3 40.14 cM
Gene Summary:	Cross-linked envelope protein of keratinocytes.[UniProtKB/Swiss-Prot Function] Transcript Variant: This variant (1) represents the longer transcript. Both variants 1 and 2 encode the same protein.