

## Product datasheet for MC209542

### Tnfrsf18 (NM\_009400) Mouse Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Tnfrsf18 (NM_009400) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Tnfrsf18
Synonyms:	AITR; Gitr
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>MC209542 representing NM_009400 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGGGGGCATGGCCATGCTGTATGGAGTCTCGATGCTCTGTGTGCTGGACCTAGGTGAGCCGAGTGTAG  
TTGAGGAGCCTGGCTGTGGCCCTGGCAAGGTTGAGAACGGAAGTGGCAACAACACTCGTCTGCAGCCT  
GTATGCTCCAGGCAAGGAGGACTGTCCAAAAGAAAGGTGCATATGTGTACACCTGAGTACCACTGTGGA  
GACCCCTCAGTGAAGATCTGCAAGCACTACCCCTGCCAACAGGCCAGAGGGTGGAGTCTCAAGGGGATA  
TTGTGTTTGGCTTCCGGTGTGTTGCCTGTGCCATGGGCACCTTCTCCGCAGGTGCTGACGGTCACTGCAG  
ACTTTGGACCAACTGTTCTCAGTTTGGATTTCTCACCATGTTCCCTGGGAACAAGACCCACAATGCTGTG  
TGCATCCCGGAGCCACTGCCCACTGAGCAATACGGCCATTTGACTGTCATCTTCCCTGGTCATGGCTGCAT  
GCATTTTCTTCTTAACCACAGTCCAGCTCGGCCTGCACATATGGCAGCTGAGGAGGCAACACATGTGTCC  
TCGAGAGACCCAGCCATTCGCGGAGGTGCAGTTGTGAGCTGAGGATGCTTGCAGCTTCCAGTTCCCTGAG  
GAGGAACGCGGGGAGCAGACAGAAGAAAAGTGTGCATCTGGGGGTGCGGTGGCCA**TGA**

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	Sgfl-MluI
ACCN:	NM_009400
Insert Size:	687 bp



[View online »](#)

**OTI Disclaimer:** Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at [custsupport@origene.com](mailto:custsupport@origene.com) or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_009400.2](#), [NP\\_033426.1](#)

**RefSeq Size:** 1007 bp

**RefSeq ORF:** 687 bp

**Locus ID:** 21936

**UniProt ID:** [O35714](#)

**Cytogenetics:** 4 E2

**Gene Summary:** Receptor for TNFSF18. Seems to be involved in interactions between activated T-lymphocytes and endothelial cells and in the regulation of T-cell receptor-mediated cell death. Mediated NF-kappa-B activation via the TRAF2/NIK pathway (By similarity).[UniProtKB/Swiss-Prot Function]