

Product datasheet for MC209553

Cd70 (NM_011617) Mouse Untagged Clone

Product data:

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product Type:	Expression Plasmids
Product Name:	Cd70 (NM_011617) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Cd70
Synonyms:	Cd27l; CD27LG; Tnfsf7; Tnlg8a
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>MC209553 representing NM_011617 Red=Cloning site Blue=ORF
	TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC GCC <mark>GCGATCGC</mark> C
	ATGCCGGAGGAAGGTCGCCCTTGCCCCTGGGTTCGCTGGAGCGGGACCGCGTTCCAGCGCCAATGGCCAT GGCTGCTGCTGGTGGTGTTTATTACTGTGTTTTGCTGTTGGTTTCATTGTAGCGGACTACTCAGTAAGCA GCAACAGAGGCTGCTGGAGCACCCTGAGCCGCACACAGCTGAGTTACAGCTGAATCTCACAGTTCCTCGG AAGGACCCCACACTGCGCTGGGGAGCAGGCCCAGCCTTGGGAAGGTCCTTCACACACGGACCAGAGCTGG AGGAGGGCCATCTGCGTATCCATCAAGATGGCCTCTACAGGCTGCATATCCAGGTGACACTGGCCAACTG CTCTTCCCCAGGCAGCACCCTGCAGCACAGGGCCACCCTGGCTGTGGGCATCTGCTCCCCCGCTGCGCAC GGCATCAGCTTGCGTGGGCGCTTTGGACAGGGCCACCTGGCCATTACAGCGCCTGACATACCTGG TCCACGGAGATGTCCTCTGTACCAACCTCACCCTGCCTCTGCTGCCGCCGCACGCGCTGACATACCTGG TCCACGGAGATGTCCTCTGTACCAACCTCACCCTGCCTCTGCTGCCGCCGCACGCTGACATGCCCTT CTTTGGAGTTCAGTGGATATGCCCTTGA
	ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT ACAAGGATGACGACGATAAGGTTTAA
Restriction Sites:	Sgfl-Mlul
ACCN:	NM_011617
Insert Size:	588 bp



This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2024 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US

Cd70 (NM_011617) Mouse Untagged Clone – MC209553

OTI Disclaimer:	Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at <u>custsupport@origene.com</u> or by calling 301.340.3188 option 3 for pricing and delivery. The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <u>More info</u>
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	 Centrifuge at 5,000xg for 5min. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. Close the tube and incubate for 10 minutes at room temperature. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>BC125407, AAI25408</u>
RefSeq Size:	744 bp
RefSeq ORF:	588 bp
Locus ID:	21948
UniProt ID:	<u>055237</u>
Cytogenetics:	17 29.69 cM
Gene Summary:	Cytokine that binds to CD27. Plays a role in T-cell activation. Induces the proliferation of costimulated T-cells and enhances the generation of cytolytic T-cells.[UniProtKB/Swiss-Prot Function]