

Product datasheet for **MC209660**

Gm4836 (NM_009529) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids

Tag: Tag Free

Symbol: Gm4836

Synonyms: 22526

Mammalian Cell Selection: Neomycin

Vector: pCMV6-Entry (PS100001)

E. coli Selection: Kanamycin (25 ug/mL)

Fully Sequenced ORF: >MC209660 representing NM_009529
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGTCTATTAAGAACTGTGGGTGATACCAAAGGATGGTTACTTATTACTTCTTGACTATGATTCAGATG
AAGAAGAAGAGCAGGCTCATTCTGAAGTAAAGAGGCCAGCATTGGCAAACATGAGAATATGCCGCCTCA
CGTAGAAGCAGATGAAGATATAAGAGATGAACAAGACAGTATGTTGGATAAATCTGGAGAAAACGTAAGT
TTCTCAGAGGAATGGCAGCGTTTTGCACGTTCTGTAGAGACCAATGGAAAACGGAATTTGTTAAGTG
GTGAACAGCAAGTTAGGAATGCTTCAGAATTGGACCTTATGGAAGTACAGAATCCAGTAACATCATGATGA
TGGGAATGCAATCCTGAAGAAGTAGTTGGAGATACAGAAAGAAGATCAACAACAAATTGTGTGAGCAG
AAGTTTGATATGGATATACAGAAATCAATGAAGAACAAGAAAAATCAGTGAACAATTATCAAAAAGAAC
AACAAGCATTAAACTGTCCGAATGTAGTCAGAGTCCGACCATGGAAGCAATTGAAGACATGCATGAGAA
GTCCATGGAGGGTTTGATGAACATGGAGACCAACAACACTACGATATGCTTTTGTAGATGGTGAAGAG
ACATTATGA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI

ACCN: NM_009529

Insert Size: 639 bp



OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_009529.3, NP_033555.1</u>
RefSeq Size:	917 bp
RefSeq ORF:	639 bp
Locus ID:	22526
Cytogenetics:	X