

Product datasheet for **MC209731**

Psma6 (NM_011968) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Psma6 (NM_011968) Mouse Untagged Clone
Tag: Tag Free
Symbol: Psma6
Synonyms: IOTA
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >MC209731 representing NM_011968
Red=Cloning site **Blue**=ORF **Orange**=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGTCCCGGGTCCAGCGCCGGTTTTGACCGTCACATTACTATCTTCTCTCCGAGGGCCGACTCTACC
AAGTAGAATATGCTTTAAAGCTATTAACCAGGGTGGACTTACATCTGTTGCTGTCAGAGGAAAGGACTG
TGCAGTTATTGTCACACAGAAGAAAGTACCCGACAACTACTGGATTCCAGCACAGTGACTCACTTGTTCC
AAGATAACGGAAAGCATTGGCTGTGTGATGACAGGAATGACAGCTGACAGCAGATCCCAGGTACAGAGGG
CACGCTATGAGGCAGCTAATTGGAATACAAGTATGGCTATGAAATTCCTGTGGACATGCTGTGCAAAAAG
AATTGCTGATTTCTCAAGTCTACACACAGAACGCTGAAATGAGGCCACTCGGTTGTTGTATGATTTTA
ATCGGTATCGATGAAGAGCAAGGCCCTCAAGTGTACAAGTGTGATCCTGCAGGTTACTACTGTGGCTTTA
AAGCCACGGCGGCAGGAGTTAAACAAACAGAGTCAACCAGCTTCCTTGAAAAGAAAGTGAAGAAGAAATT
TGATTGGACATTTGAACAGACAGTGAAACTGCAATTACATGCCTGTCTACTGTTCTGTCGATTGACTTC
AAACCTTCAGAAATAGAAGTTGGAGTAGTTACAGTTGAAAATCCTAAATTCAGGATCTTACTGAAGCAG
AGATTGACGCTCACCTTGTGCTCTAGCAGAGAGAGACT**GTA**

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_011968
Insert Size: 741 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_011968.3 , NP_036098.1
RefSeq Size:	997 bp
RefSeq ORF:	741 bp
Locus ID:	26443
UniProt ID:	Q9QUM9
Cytogenetics:	12 C1
Gene Summary:	<p>Component of the 20S core proteasome complex involved in the proteolytic degradation of most intracellular proteins. This complex plays numerous essential roles within the cell by associating with different regulatory particles. Associated with two 19S regulatory particles, forms the 26S proteasome and thus participates in the ATP-dependent degradation of ubiquitinated proteins. The 26S proteasome plays a key role in the maintenance of protein homeostasis by removing misfolded or damaged proteins that could impair cellular functions, and by removing proteins whose functions are no longer required. Associated with the PA200 or PA28, the 20S proteasome mediates ubiquitin-independent protein degradation. This type of proteolysis is required in several pathways including spermatogenesis (20S-PA200 complex) or generation of a subset of MHC class I-presented antigenic peptides (20S-PA28 complex).[UniProtKB/Swiss-Prot Function]</p> <p>Transcript Variant: This variant (1) represents the longer transcript and encodes the longer isoform (1).</p>