

Product datasheet for MC209757

Spo11 (NM_012046) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Spo11 (NM_012046) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Spo11
Synonyms:	AI449549; Spo11a; Spo11b
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>MC209757 representing NM_012046 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGCGTTCGCGCCTATGGGGCCGAGGCCTCGTTCCTCGACGCCCTGGATCGGCACAGGGCTTCCTGT
TGGCCATGGTGAAGAGAGGGCGCAGGGGAGACCCCTGCCGGGGCCACCCCGTGGCCTCTAGTTCGAGGT
TCTTACAGCTATAGAAAATATTCCAAGACATAATCAAAGCTTGGCAAGAAATGAAGTGCCTGCCTTC
ACAATAGACAACAGATCGAGCTGGGAAAACATAATGTTTGATGATTCTGTGGCCTTCGGATGATACCTC
AGTGTACCACAAGAAAAATCAGAAGCGATTACCAAAAATCAGTTAAGAAATTCGCTCTGATTCTGAAAGT
ATTGTCCATGATTATAAATTAATACAGAGCGACACTTATGCAACCAAGAGAGACATATACTACACTGAC
AGCCAGCTCTTTGGGAACCAGGCTGCGGTGGACAGCGCCATCGATGACATTTCTGTATGCTGAAAGTGC
CCAGGAGGAGTCTGCACGTGCTATCTACTTCCAAGGGATTGATTGCTGGCAACTTGAGATACATGGAGGA
AGATGGTACCAGAGTCCAGTGTACCTGTAGTGCCACGGCTACTGCTGTGCCGACTAACATTCAAGGAATG
CAGCATCTGATCACAGATGCGAAGTTTCTGTTAATAGTCGAGAAGGATGCAACATTTTCAGCGGCTCCTGG
ACGACAACCTCTGCAGCAGGATGTCCCGTGCATCATGGTTACGGGAAAGGGCGTTCAGATCTGAACAC
GAGGCTCTTGGTAAAGAAGCTGTGGGACACCTTCATATTCCTGTTTTACACTGGTTCGATGCAGATCCC
TACGGCATCGAGATAATGTGCATCTATAAGTACGGATCCATGTCCATGTCTTTTGAAGCTCACAATCTCA
CTATCCCAACAATCAGATGGCTTGGTCTCCTCCCTTCTGATATCCAGAGGTTAAATATACCTAAGGATAG
TTTGATCCCACTGACAAAGCATGACCAGATGAAGCTGGACAGCATCCTGAAGAGGCTTACATTACCTAC
CAGCCCTCTGGAAAAAGAGCTGGAGATGATGGCAGACTCTAAGATGAAGGCAGAGATCCAAGCCTTGA
CCTTGCTGCTGTCAGACTACCTGTCCAGAGTGTACTTACCAACAAGCTGAGTTTGGAGGATGGATCTA
A

ACGGTACGGCGGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online >](#)

Chromatograms:	https://cdn.origene.com/chromatograms/ja2262_c02.zip
Restriction Sites:	Sgfl-Mlul
ACCN:	NM_012046
Insert Size:	1191 bp
OTI Disclaimer:	<p>Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.</p> <p>The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info</p>
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_012046.2 , NP_036176.1
RefSeq Size:	1716 bp
RefSeq ORF:	1191 bp
Locus ID:	26972
UniProt ID:	Q9WTK8
Cytogenetics:	2 95.64 cM
Gene Summary:	<p>Isoform 1: Component of a topoisomerase 6 complex specifically required for meiotic recombination. Together with TOP6BL, mediates DNA cleavage that forms the double-strand breaks (DSB) that initiate meiotic recombination (PubMed:26917764). The complex promotes relaxation of negative and positive supercoiled DNA and DNA decatenation through cleavage and ligation cycles. Essential for the phosphorylation of SMC3, HORMAD1 and HORMAD2 (PubMed:22346761).[UniProtKB/Swiss-Prot Function]</p> <p>Transcript Variant: This variant (1) represents the longest transcript and encodes the longest isoform (a).</p>