

Product datasheet for MC209764

Tspan32 (NM_020286) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Tspan32 (NM_020286) Mouse Untagged Clone
Tag: Tag Free
Symbol: Tspan32
Synonyms: Art-1; AW208513; BB235973; D7Wsu37e; Phemx; Tssc6
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Fully Sequenced ORF: >MC209764 representing NM_020286
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGGCCACCATGGTGGTTGTCATTCACCTTTGGAGACCACTTTACTGTGATTGGCCATGCATCGTTGGAGA
 GGAACCCCTATGAAACACTGCGCTATTGGGCCTTCTATGTGGGGATCAGCCTAGCGGGGCTCCTGAGCCT
 GGGAGCTGCTCTGAGCACCATAGCCACTGTGAGGGAGGCCATGGCCTCATGGCTGCGGGTTTCCTGTGC
 TTTGCCCTGTATTCTGCATCCTGGTGCAGGTGGCCTTCTGGAGATTCTACAACCCACCCAGGTGGAGG
 ATGCGGTGTTGGACACCTATGACTTCGTGTATGATCAGGCAATGAAGAGCCCATCTAGCAACTGGTGGCA
 GGAGCTGGCCGTCATCCAGGACACGTTTCTGTGTTGTGGGAAGAAATCTCCTTTTGGGCTTCTGGTGAGC
 ACTGGAGCCATCATGTGTCAGGGCCGGGAGGCCATGAGAGAGGACTGCCTCCAGAGCATCAGGAACGTTT
 TGTGGACACACTACAGTATTGCTTCCATCCTGACCTGCACAGCCTGGCCCTCACGGTATATGCTATGAT
 GCTCTGTGCCTTCTGTGGTTTGCATTCACTCTACCATGGCTTGGATAGAAAAGGCAGATACTCCCTG
 ACACCACCAAGATCCCATGGTTTCCAGACCCAGGAACCCAGCCTCTTCAGATGGAC**TAG**

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_020286
Insert Size: 690 bp


[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	NM_020286.3 , NP_064682.1
RefSeq Size:	1697 bp
RefSeq ORF:	690 bp
Locus ID:	27027
UniProt ID:	Q9JHH2
Cytogenetics:	7 88.09 cM