

Product datasheet for **MC209793**

Skp2 (NM_013787) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Skp2 (NM_013787) Mouse Untagged Clone
Tag: Tag Free
Symbol: Skp2
Synonyms: 4930500A04Rik; FBXL1; FWD1; p45
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >MC209793 representing NM_013787
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**GCGATCGCC**

ATGCATAGGAAGCACCTTCAGGAGATCCGGACCAGAGTGCCAACGTACCACCAGCTTCACGTGGGGAT
 GGGATTCCAGCAAGACTTCTGAAGTCTATCAGGCATGGGTGTCTCGGCCTTGGAGAAGGAGGAGGTGGA
 CAGTGAGAACATCCACATGGACTGCTCTCAAACCTCGGCCACCCAGAGCCCTCCAAGAAACGAGTC
 AAGGGCAAAGGGAGTGACAAAGACTTTGTGATCATCCGTCCGGCCGAAGCTTAGTCGGGAGAACTTTCCAG
 GTGTCTCCTGGGACTCCCTCCAGATGAGCTGCTCCTTGGGATCTTTTCTGTCTGTGCCTCCCTGAGCT
 GCTGAGAGTCTCGGGCGTTTGAAGAGGTGGTACCGCCTCTCGCTCGATGAGTCTCTTGGCAGTCCCTC
 GACCTCGCGGGTAAAAATCTGCACCCAGACGTGACTGTGCGCTTGCTCTCCCGGGGGTGGTGCCTTCC
 GCTGCCCTCGGTCCTTTATGGAGCAGCCGCTGGGTGAAAGCTTCAGCTCTTCCGGGTACAGCACATGGA
 CCTGTGAACTCAGTGATAAATGTGTCGAACCTCCATAAGATTCTGTCCGAGTGTCCAAGCTGCAGAAT
 CTAAGCCTGGAAGGCCTGCAGCTCTCAGACCCATTGTCAAGACTTTGCACAGAATGAAAACCTTGGTGC
 GACTAAACCTTTGTGGTGCTCTGGGTTTCTGAATCTGCCGTGGCGACTCTGCTAAGCAGCTGCTCCAG
 ACTGGATGAGCTAAATCTCCTGGTGCTTTGACTTCACTGAAAAGCACGTGCAAGCGGCTGTGGCACAT
 TTACCAAACACCATCACCCAGCTGAACCTCAGCGGCTACCGAAAAGAACCTCCAGAAAACAGATCTTTGTA
 CCATAATTAACGATGCCCAACCTCATCCGCTCGACTTAAGTGACAGTATCATGCTAAAGAATGACTG
 CTTTCCAGAATTTTTCAACTCAACTACCTCCAACACCTCTCGCTCAGCCGGTGTATGATATAATACCT
 GATACTCTACTTGAAGTGGAGAAATTCCTACACTAAAAACGCTACAAGTTTTTGGAAATCGTGCCAGAGG
 GAACCTTCAGCTACTGAGGGAAGCTTCTCCTCGGCTGCAGATTAAGTGCCTATTTACCACCATTGC
 AAGGCCAACTATGGACAGCAAGAAGAACCTGGAGATTTGGGTATCAAGTGCCGACTGACTCTGCAAAG
 CCCAGTTGTCTA**TGA**

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Restriction Sites:	Sgfl-Mlul
ACCN:	NM_013787
Insert Size:	1275 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_013787.3 , NP_038815.1
RefSeq Size:	3210 bp
RefSeq ORF:	1275 bp
Locus ID:	27401
UniProt ID:	Q9Z0Z3
Cytogenetics:	15 A1
Gene Summary:	<p>Substrate recognition component of the SCF (SKP1-CUL1-F-box protein) E3 ubiquitin-protein ligase complex which mediates the ubiquitination and subsequent proteasomal degradation of target proteins involved in cell cycle progression, signal transduction and transcription. The SCF complex provides substrate specificity and interacts with both, the E2 ubiquitin-conjugating enzyme and the substrate. Specifically recognizes phosphorylated CDKN1B/p27kip and is involved in regulation of G1/S transition. Degradation of CDKN1B/p27kip also requires CKS1. Promotes ubiquitination and destruction of CDH1 in a CK1-Dependent Manner, thereby regulating cell migration (By similarity).[UniProtKB/Swiss-Prot Function]</p> <p>Transcript Variant: This variant (1) represents the longer transcript and encodes the longer isoform (a).</p>