

## Product datasheet for MC209799

### Klra15 (NM\_013793) Mouse Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Klra15 (NM_013793) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Klra15
Synonyms:	Ly490
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>MC209799 representing NM_013793 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGGATCGCC**

ATGAGTGAGCAGGAGGTCACTTATTCAATGGTGAGATTCATAAATCTGCAGGATTGCAGAAACAAGTGA  
GACCTGAGGAGACTAAAGGGCCAGAGAAGCTGGCTACAGAAGGTGTTCCACTGGAAGTTCATTGT  
GATAGCTCTTGGCATCTTCTGTTTCTTCTGGTAGCTGTTTCAGTGTGGCAATAAAAAATTTTCAG  
TATGATCAACAAAAAACTGCAGGAAATCTAAACCACCACAATAACTGCAGCAACATGCAAAGTGACA  
TCAACTTGAAGGATGAACTGCTGAAAAATAAGTCCATAGAGTGTGATCTTCTGGAATCCCTCAACAGGGA  
TCAGAACAGATTGTATAATAAAACCAAGACTGTTTTAGATTCCTTAAAGCACACAGGCAGAGGTGTTAAA  
GTATACTGGTTCTGCTATGGTATGAAATGTTATTTTCGTCATGGACAGAAAACCATGGAGTAGATGTA  
AACAGAGCTGCCAGAATTCAGCTTATTCTTCTGAAGATAGACGATGAGGATGAACTGAAGTTCCTTCA  
GCTCGTGGTTCCTTCAGACAGTTGCTGGATTGGATTGTCATATGATAATAAGAAAAAAGATTGGGCATGG  
ATTGAAAATGGCCCATCTAACTTGCCCTGAACACAACGAAATACAATATAAGAGATGGAGGATGTATGT  
TCTTATCTAAAAGTACTAGACAATAATTACTGTGATCAATCATTATCTGTATTTGTGGGAAGAGATT  
GGACAAATCCCTCATTGA

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	Sgfl-MluI
ACCN:	NM_013793
Insert Size:	789 bp



[View online »](#)

<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">NM_013793.2</a></u> , <u><a href="#">NP_038821.2</a></u>
<b>RefSeq Size:</b>	1083 bp
<b>RefSeq ORF:</b>	789 bp
<b>Locus ID:</b>	27423
<b>Cytogenetics:</b>	6 F3