

Product datasheet for **MC209817**

Hspb7 (NM_013868) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Hspb7 (NM_013868) Mouse Untagged Clone
Tag: Tag Free
Symbol: Hspb7
Synonyms: 27kDa; cvHsp; Hsp25-2
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >MC209817 representing NM_013868
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGAGCCATCGGACCTCCTCCGCCTTCAGAGCGGAGAGAAGCTTCGGTTCCTCTCCTCCTCCTCCTCCT
CCTCTTCTTCTCGGCCCTCCCGTGCCCTCCAGCCAGGACCCACCCATGGAGAAGGCTTTGAGCATGTT
TTCAGACGACTTTGGCAGTTTTATGCTGCCCACTCGGAGCCCTGGCCTCCCAGCCCGTCCCGGGGGC
CAAGGAAACATCAAGACCCTCGGGGATGCCTACGAGTTTACAGTGGACATGAGAGACTTCTCCCCTGAAG
ACATCATTGTCAACCCTTCAACAACCACATCGAGGTGCGGGCTGAGAAGCTGGCAGCTGATGGCACAGT
CATGAACACCTTTGCTCACAAGTGCCAACCTACCAGAGGATGTGGACCAACATCGGTGACCTCGGCCCTT
CGAGAGGATGGTAGCCTCACCATCCGGGCGCGGCCACCCCCACACAGAACATGTCCAGCAGACCTTCC
GGACTGAGATAAAAAT**TGA**

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_013868
Insert Size: 510 bp
OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_013868.4</u> , <u>NP_038896.2</u>
RefSeq Size:	2768 bp
RefSeq ORF:	510 bp
Locus ID:	29818
UniProt ID:	<u>P35385</u>
Cytogenetics:	4 D3