

## Product datasheet for **MC209846**

### Siva1 (NM\_001161737) Mouse Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Siva1 (NM\_001161737) Mouse Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Siva1  
**Synonyms:** CD27bp; Siva  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**Fully Sequenced ORF:** >MC209846 representing NM\_001161737  
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGGATCGCC**

ATGCCAAGCGGAGCTGCCCGTTCGCAGACGCAGCCCCGCTCCAAGTCCACGTGGGCTGAAAG  
AGCTGAGCCACGGTGTGTTCCGCGAGCGCTACTCAGCGAGGTCTTCGGCCTTCCAGGACAGCGCCAT  
CGCTTGTTCATCGTGCATGAGATCTGTGGATGGGAAGCGGTCTGCAGCCAGTGCAGCGGGCCCTGTGT  
GGCAGTGTGTATACACCTGCTGGGGCTGCGGTGCTTTGGCCTGTGTGCTGTGTGGCCTTGCAGACTATG  
CCGACGATGGTGAGAAGACACTGTGCACCAGCTGTGCTATGTTTGAAGCCTGA

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** SgfI-MluI

**ACCN:** NM\_001161737

**Insert Size:** 333 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).



[View online >](#)

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001161737.1](#), [NP\\_001155209.1](#)

**RefSeq Size:** 587 bp

**RefSeq ORF:** 333 bp

**Locus ID:** 30954

**UniProt ID:** [O54926](#)

**Cytogenetics:** 12 F1

**Gene Summary:** Induces CD27-mediated apoptosis. Inhibits BCL2L1 isoform Bcl-x(L) anti-apoptotic activity. Inhibits activation of NF-kappa-B and promotes T-cell receptor-mediated apoptosis (By similarity).[UniProtKB/Swiss-Prot Function]  
Transcript Variant: This variant (2) lacks an exon in the coding region, compared to variant 1. The encoded isoform (2) is shorter, compared to isoform 1.