

Product datasheet for MC209901

Ifi27 (NM_194066) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Ifi27 (NM_194066) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Ifi27
Synonyms:	1110013J02Rik; 2900026P10Rik; D12ErtD647e; Ifi27I1; ISG12a
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin
Fully Sequenced ORF:	>MC209901 representing NM_194066 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGCGTTTTCGGGCACAGGGACACTGGTGGCTTCCATTGTCTCCAAGATGACATCTTCAGCAGCCATGG
TCAAAGCTGGAGGGCAGTCTCTGGCACCATTCTAGCTGGTTCGAAGGGGCTCATCTGTTAACACAGAG
TGCCCTGGGCTCTGCAACATCTGCACTTGGAGCACTGAAGGTTGGCACCATTTATCTGGGTTCTCTGCC
TCCACCTTGGCTGCTTCCCCATTGGAGAGGAACCTGCTTTGGTGGCAGGGAGAAGAAAAGTGACTTCA
GATTTTCTTCCCGCAGCAGGCAAGCAGCAAGAAGGCAGAGAAAGGAAGGACAAAGCAAAATCTACCT
GAAGGGTTTTGAGTTCTCCAGAGCCATGACAGTAGCTGCTGTACCACCTGCGCTGAGTGCTGTGGGCTTC
ACCGCCTCAGGAATTGCAGCCTCTCTAGCAGCTAAGATGATGTCCTTGTGAGCTATTGCTAATGGGG
GTGGAGTCCCAGCTGGTGGCCTGGTGGCCATTCTGAGTCTGCTGGAGCTGCTGGACTCTCCGTGCCATC
TACTGTATCGTGGGCTCTGCAGGCTGCGGTTGTGGCCAGTGTGATGAACATCTGTGAGAGTTTCTAC
CCTTTTCTCATGGGATCAGAGGTGGCAGACATGGCGACAGAGGTGGCAGATATGGCGATAGAGGTGGCAG
ACATTTCAACAGAAGAGGCACTGCCCACTTGTCTAACACAGAGAAGGATTAA

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	Sgfl-MluI
ACCN:	NM_194066
Insert Size:	753 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_194066.2 , NP_919042.2
RefSeq Size:	1023 bp
RefSeq ORF:	753 bp
Locus ID:	52668
Cytogenetics:	12 52.93 cM