

Product datasheet for **MC209909**

Cd300lg (NM_001160712) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Cd300lg (NM_001160712) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Cd300lg
Synonyms:	2310016B05Rik; Clm9; D11ErtD736e
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>NM_001160712.1 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGAGGCCTCTGGTCTGCTATGGGGCTGCCTGGTGTCCAGGTTATGAAGCCCTGAAGGGTCCAAAGG
AGATCAGTGGATTTGAAGGTGACACCGTGTCCCTGCGGTGTACCTACGTGGAGAAGATGAAGGAGCACAG
GAAGTATTGGTGCCGGCAGGGTGGCATCCTGGTGTACGCTGCGGTGACATTGTCTACGCAATCAGGAC
CAGGAGGTGACTCGAGGCAGGATGTCCATCCGAGACAGTCCCAAGAGCTCTCGATGACCGTGATCATGA
GGGACCTTACCCTGAAGGATTACAGGAAGTACTGGTGTGGGATTGACAGACTGGGCCGCGATGAGTCTTT
TGAGGTTACACTATTGTCTTTCCAGGAAAAGGGACAATGCTGTCCCTGCTGGCACCTGCTGTCCTTCT
TCCCCACCCCTCCTTCCAGCCTTTACACCTACAAGGAGCCTGCAGCCCAAGGCAAAAGCCTGGCAGA
CTCAGCTCCCAGAACCGAGGAGCTCCCCTCCAGTCTGCTGGTGTCCCTTACCACGCCACAGGACTCCAG
GGCTGTAGCCAGCAGTGTCTCCAAGCCAGTGTCCATCCCGATGGTCCGCATGATGGCCCCAGTCTGT
ATACTCTGTCCCTGCTGTTGGCTGCAGGACTAATTGCCTTTGGCAGCCACATGCTCCGGTGGAGAAAGA
AAGCTTGGCTGGCCACAGAGACACAGAAGAACGAGAAGGTCTACCTTGAACCTCGCTGCCAGGGAACGG
CTGGACCACTGAAGACTCGACGATAGACCTTGCAGTGACTCCTGAATGTCTCAGAAACCTCAACCTTCT
GCTGTGCCCTCTCCTGAGACACAGAACCCTCAGTCAGTCTACAGAGGAGGAAGAGGAGCTCGTTCCTGG
ACGACGACAAGGAGGACGTGATGGCACCCCTCCCTTGCAGATGTCTGCCGAGGAAGTGGCCCTTCTCTGA
GTTTCATCTCTGTGTA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCTGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Chromatograms: https://cdn.origene.com/chromatograms/ja1435_a10.zip



[View online »](#)

Restriction Sites:	Sgfl-Mlul
ACCN:	NM_001160712
Insert Size:	996 bp
OTI Disclaimer:	<p>Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.</p> <p>The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info</p>
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001160712.1 , NP_001154184.1
RefSeq Size:	2480 bp
RefSeq ORF:	996 bp
Locus ID:	52685
UniProt ID:	Q1ERP8
Cytogenetics:	11 65.48 cM
Gene Summary:	<p>Receptor which may mediate L-selectin-dependent lymphocyte rollings. Binds SELL in a calcium dependent manner. Binds lymphocyte.[UniProtKB/Swiss-Prot Function]</p> <p>Transcript Variant: This variant (3) lacks an in-frame segment, as compared to variant 1. The resulting isoform (3) lacks an internal segment, as compared to isoform 1.</p>