

Product datasheet for MC210106

Tssc4 (NM_001115085) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Tssc4 (NM_001115085) Mouse Untagged Clone
Tag: Tag Free
Symbol: Tssc4
Synonyms: AA241958; ESTM671070
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >MC210106 representing NM_001115085
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGGCGGAGACAGAGGCTGGCCTGGAGGCTGATGAACCCACCGAGGATGACACCCTGCCTTCTGACACCG
TCTCCCTCAGCGACTCGGACTCGGACCTTAGCTTACCTAGTGGTGGAGGTCCAAGTGTGTCCCAGAG
GAGGCTTTCGGAGGAGGGCCAGGAGGACTCCGGCCCTGAGGACCCTCCATCGCCCCACGGGCACCCCTC
ACCACTGCAGTGCAGCCATTCCACCTCCGAGGCATGAGTTCTACCTTCTCCAGCGCAGCCACAGCATCT
TTGATTGTCTGGAGAGTGCAGCCCGGCAGGCACCGTGTCTGCTCCCAAACAGTGTGAGTGACAATGG
CAGCTTACGGCGGCCCGTGACTCCCCAAGTCAGACTCCAGCAAGGGGCCTGAGCAGAGTGCATGGGAAC
ACTGGCCCGACTAGGGTACTCCAGTCCCTGACTATGTGTACACCCTGAGCGTTGGACCAATACAGTC
TGGAGGATGTGTCTGAAGCCAGTGAAGCAGCAATCGTGATGCTGCCCTGGCCTTCTGAGCTCTCGAAG
CCAGGTGTCCACACAGACTATGTGCCTTCTTTCAACCAGGACCCTCTAGCTGTGGGGAGGGAAGAGTG
GTCTTACCAAACAGTGAAGGACAGTGAAGGCTAGGGCTGAGAGGAAGAGGGTCTAAAGAAGGGGGTTG
GGTCAGGTGCGGGGGTGAAGCCGCTGTGAGCTAGCCATCTCGTGGACCTGAGGCTGAGGAATGGAG
TGGCCACAGGGACAGCCAGAGGTGGTGTACCTCAGAAGCTGCCATCCTGAGTCTCATCAGGCCCC
ATAGGTGTGAAGGCTGTTGGTTTCCATGGCAGCAAGAAGCGGAGCAGAGACCCTCCGGAACAGGGACG
GTAACCCTGGGGGCGGGGTCTGAGAGAGGACCCTCAGTT**TGA**

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_001115085
Insert Size: 954 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001115085.1</u> , <u>NP_001108557.1</u>
RefSeq Size:	1364 bp
RefSeq ORF:	954 bp
Locus ID:	56844
UniProt ID:	<u>Q9JHE7</u>
Cytogenetics:	7 88.11 cM