

Product datasheet for MC210207

Mrpl38 (NM_024177) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Tag:	Tag Free
Symbol:	Mrpl38
Synonyms:	1110036N21Rik; 4733401F03Rik; MRP-L3; R75240; Rpm13
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>MC210207 representing NM_024177 Red=Cloning site Blue=ORF Orange=Stop codon

TTTGTGAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGGCGGCGCCCTGGTGGCGTGC GGCGT TTTTCGGAATTGGGAGATGCCGGGGCTTCAGCACCTCGGCCT
TCCTCAGCCGCCGGACACCCCGCTCGGGCCCATGCCCAACGAGGACATCGACGTGAGCAACCTGGAGCG
GCTGGAGAAGTACCGGAGCTTCGAGCGGTACCGCCCGGGCGGAGCAGGAGGCGCGCCCGCACTGG
TGGCGGACCTACCGGGAGCATTTCTGTAAGGAGACAGATCCCAAAGACAAAATTGACATTGGGCTACCCC
CACCCAGAGTCAGTCGGACCCAAAAGCTACTGGAGCGGAAGCACTTTCTGCGGGAACCTTCGTGCCAATGT
GGAAGAAGAACGGGCAGCCCGCTCCGCACGGCTAGCATCCCACTGGAGGCAGTGCGGGCAGAGTGGGAG
AGGACGTGTGGCCCTACCACAAGCAGCGTCTGGCTGAGTACTATGGCCTCTACCGAGACCTGTTCCACG
GTGCCACTTTTGTGCCCTGGGTCCCTCTACACGTGGCTTACGCCGTGGGAGAAGAGGACCTGATCCCGT
GTATCATGGCAATGAAGTGACTCCGACCGAGGCTTCCCGGGCCCCAGAGGTGACCTATGAGGCAGACAAG
GATTCCTTGTGGACACTGCTATTCACTTGGACGGACACCTGCTGGAGCCGGACGCTGAGTATGTTT
ACTGGCTGCTAACCAACATCCCAAGCAACAGAGTGGCTGAAGGACAGGAGACATGTCCTACCTCCCTCC
TTTCCCTGCTCAGGCTCTGGCTTCCACCGCTTTGCCCTTCTTGCTCTTCAAGCAGGACAAGCCAATCAAC
TTCTCAGAGGACACCCGTCCTCACCCTGCTACCAGCTGGCCAGAGGACCTTCCGCACTTTTGATTCT
ACAAGAGACACCAGGAAGCCATGACTCCAGCCGGTCTGGCTTTTTTCCAGTGCCGCTGGGATGACTCTGT
CACCCACACCTTCCACCAGCTTCTAGACATGCGGGAACCTGTGTTGAATTTGTACGGCCACCCCTTAC
CACCCCAAACAGAAGCGCTTCCCCACGAGCAGCCCTACGCTATCTGGACCGATACAGGGACAGTCATG
AGCCACCTATGGAATCTACTGA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



Restriction Sites:	Sgfl-Mlul
ACCN:	NM_024177
Insert Size:	1143 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_024177.3, NP_077139.2</u>
RefSeq Size:	1411 bp
RefSeq ORF:	1143 bp
Locus ID:	60441
UniProt ID:	<u>Q8K2M0</u>
Cytogenetics:	11 E2