

## Product datasheet for **MC210278**

### Josd2 (NM\_025368) Mouse Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Josd2 (NM\_025368) Mouse Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Josd2  
**Synonyms:** 1110007C05Rik  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**Fully Sequenced ORF:** >MC210278 representing NM\_025368  
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**GCGATCGCC**

ATGTCTCAGGCTCCGGAAGCAAGCCGAGTCCACCCTCAGTGTATCATGAGCGGCAGCGGCTAGAGCTGT  
GTGCTGTCCATGCTCTCAACAACGTCTTCAAGAGCAGCTCTTTAGCCAGGAGGCTGCCACGAAATCTG  
CAAGAGGCTAGCCCCAGATTCCCGGCTGAACCCCATCGCAGTCTTTGGGCACCGGCAACTATGATGTC  
AACGTGATCATGGCTGCCCTGCAGGGCCTGGGCCCTGGCCGCTGTGTGGTGGGACAGGAGGAGACCCCTGT  
CCCAGCTGGCCCTGCCCCAGTACTAGGCCTGATCTTGAACCTACCCTCTCCTGTGTCGCTGGGGCTGT  
GTCCCTGCCACTGCGCCGCCGACTGGGTGGCCCTGCGACAGGTAGATGGCATCTACTATAATCTGGAC  
TCAAAGCTCCGGGCCCTGAAGCCCTGGGGGACGAGGATGGAGTCAGGACCTTCTGGCAGCTGCCTGG  
CCAAGGCCTGTGCGAGGTGCTATTGGTGGTGACCAAGGAGGTGGAGGAGCCGGCTGCTGGCTGAACAC  
AAGCTGA

**ACGCGT**ACGCGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** SgfI-MluI  
**ACCN:** NM\_025368  
**Insert Size:** 567 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u>NM_025368.4, NP_079644.1</u>
<b>RefSeq Size:</b>	936 bp
<b>RefSeq ORF:</b>	567 bp
<b>Locus ID:</b>	66124
<b>UniProt ID:</b>	<u>Q9CR30</u>
<b>Cytogenetics:</b>	7 B3
<b>Gene Summary:</b>	<p>Cleaves 'Lys-63'-linked poly-ubiquitin chains, and with lesser efficiency 'Lys-48'-linked poly-ubiquitin chains (in vitro). May act as a deubiquitinating enzyme (By similarity). [UniProtKB/Swiss-Prot Function]</p> <p>Transcript Variant: This variant (3) differs in the 5' UTR, compared to variant 1. Variants 1, 2, 3, 4 and 5 all encode the same isoform (a). Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.</p>