

Product datasheet for MC210280

Znrd1 (NM_023162) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Znrd1 (NM_023162) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Znrd1
Synonyms:	1110014N07Rik; Rpa12
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>MC210280 representing NM_023162 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCCGCATCGCC

ATGGAGCTGGCCCGCCCCGCTCCAACCTTTCAATCAGACTTGGATTTCTGCCAGATTGTGGCTCAGTCC
TGCCGCTGCCTGGAATTCAGGACACTGTCATATGTTTCGCGTTGTGGCTTCTCCATCGATGTGCGAGATTG
TGAAGGGAAGTTGTGAAGACCTCAGTAGTGTTCAACAAGCTGGGGCTACCATACTTTGTCTGTGGAT
GAGGGACCTGAATTACAGGGACCACTGATTGACAGGCGCTGCCCTCGATGTGGTCACGAGGGAATGGCAT
ATCACACCAGACAGATGCGCTCAGCTGATGAAGGACAGACGGTCTTCTATACCTGTATCAACTGCAAGTT
TCAGGAAAAGGAAGATTCTGA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	Sgfl-MluI
ACCN:	NM_023162
Insert Size:	372 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).


[View online »](#)

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_023162.4, NP_075651.1</u>
RefSeq Size:	686 bp
RefSeq ORF:	372 bp
Locus ID:	66136
UniProt ID:	<u>Q791N7</u>
Cytogenetics:	17 B1
Gene Summary:	DNA-dependent RNA polymerase catalyzes the transcription of DNA into RNA using the four ribonucleoside triphosphates as substrates. Component of RNA polymerase I which synthesizes ribosomal RNA precursors.[UniProtKB/Swiss-Prot Function]