

Product datasheet for MC210303

Lce1m (NM_025420) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Lce1m (NM_025420) Mouse Untagged Clone
Tag: Tag Free
Symbol: Lce1m
Synonyms: 1110059L13Rik; AI131654; Lce5a; Sprrl10
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Fully Sequenced ORF: >MC210303 representing NM_025420
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGTCCTGTCTCAGCAGAGCCAGCAACAGTGCCAGCCCCCTCCCAAATGCAAGACCCCAAAGTGCCCTCAA
 CCTCTTCTGCTGTAGCCTGGTCTGGGGCTGCTGTGGTCCAGCTCCGGGGCTGCTGCAGCTCTGG
 ATGTGGTGGTGTCTGAGCCACCACAAGCCCCGATTCTCTCCGTCGCCGACGTACAGCTCTGGA
 TGCTGTAGCAGCGTGGCAGCAGTTGCTGTGGCAGCAGCGTGGCAGCAGTTGCTGTGGCAGCAGCGGT
 GCAGCAGCGTGGCAGCAGCGTGGCAGCAGCGTGGCAGCAGTTGCTGTGGCAGCAGCGTGGCAGCAG
 CGGTGGCAGCAGTTGCTGTGGCAGCAGCAGCGTGGCAGCAGTTGCTGTGGCAGCAGCAGTTGGTGGCAGT
 AGCGGTGGCAGTAGCGGTGGCAGCAGCTGCTGTGGTAGCAGTGGGGCAGCAGTTGGTGGCAGCAGCAGT
 GCAGCAACAGCCATCAGCAGTCTGGGGCTCTAGCTGCTGCTCTGGAGTTGCTGCT**GA**

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_025420
Insert Size: 549 bp
OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_025420.2, NP_079696.1</u>
RefSeq Size:	877 bp
RefSeq ORF:	549 bp
Locus ID:	66203
Cytogenetics:	3