

Product datasheet for MC210332

1700020L24Rik (NM_025492) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: 1700020L24Rik (NM_025492) Mouse Untagged Clone
Tag: Tag Free
Symbol: 1700020L24Rik
Synonyms: RP23-249K18.3
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Fully Sequenced ORF: >MC210332 representing NM_025492
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGATAAACATGGCGTGAAGACTCCTTTGTGGAGGAAGGAGGTGGAGGATCCCGAGGCCAGGGAGGAGG
ATCTGGAGGACGATTCATCGTCCTCGTCTTCGTCGTCGGTGGAGAGAGGTGAGACCCCGAGAGCGC
GACCGAGACGGAGGAGGACTCCCGAGACGCGGAAGAGCGCGAGGCCCGCTCTGTGTCTACAGTCCGCTG
CGCCAGGAGTCCAGCAGCCAGCAGGTGGCGCTGCTGCGACGCTCAGACAGCAGCTTCTGGGGTTGGCTTA
GCCCGTTCGCTCTGCTGGGTGGCCTGGCGGCCCGGCAGACAGGAAGCGGGTGCCTCCGGAGGAACCGTG
CGTGTGGAGACGCGGCCGCGCCACCGCGCCGGGTGGCTGTGCGCGCTGTGAGATTCTTTCTGCAAG
AAATGCAAGACTCTGCACAGCCACCCGGCCTACGTGGAGCACTGTATTCTTGAGCACCCCTGACCTGGGTA
AGGCAGAGGCCACCGCAACTCTGAGCTCATTGACTCCAGCCTCCGAGCCCGCAATGCTCCAAACTCTT
CTATCTC**TAG**

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: Sgfl-MluI
ACCN: NM_025492
Insert Size: 570 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_025492.3, NP_079768.2</u>
RefSeq Size:	931 bp
RefSeq ORF:	570 bp
Locus ID:	66330
UniProt ID:	<u>Q8C1R3</u>
Cytogenetics:	11 C