

Product datasheet for MC210354

Krtap3-3 (NM_025524) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Krtap3-3 (NM_025524) Mouse Untagged Clone
Tag: Tag Free
Symbol: Krtap3-3
Synonyms: 2310020D06Rik; KAP3.3; KRTAP3.3
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >MC210354 representing NM_025524
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCGCGATCGCC

ATGGCTTGCTGTGTTGCCGCTGCTGCAGTGTCCCACCGGCCCTGCCACCACCATCTGCTCCTCAGACA
 AGTCCTGCCGCTGTGGAGTCTGCCTGCCAGCACCTGCCACACACCATTGGCAGCTGGAGCCACCTG
 CTGTGACAACTGCCCCCACCCTGCCACATCCCTCAGCCCTGTGTGCCACCTGCTTCTGCTCAACTCC
 TGCCACCAACCCAGACCTGCTGACTGTCAACCTCACCACCTATGTGCAGCCAGGCTGTGAGGAGCCCT
 GTGTCCAAGGTGCTGTGA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: Sgfl-Mlul
ACCN: NM_025524
Insert Size: 300 bp
OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).



[View online »](#)

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_025524.2](#), [NP_079800.1](#)

RefSeq Size: 732 bp

RefSeq ORF: 300 bp

Locus ID: 66380

UniProt ID: [Q9D638](#)

Cytogenetics: 11 D