

Product datasheet for **MC210376**

Rpl15 (NM_025586) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Rpl15 (NM_025586) Mouse Untagged Clone
Tag: Tag Free
Symbol: Rpl15
Synonyms: 2510008H07Rik
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >MC210376 representing NM_025586
Red=Cloning site **Blue**=ORF **Orange**=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGGGTGCATACAAGTACATCCAGGAGCTATGGAGGAAGAAGCAGTCCGACGTGATGCGCTTTCTTCTGA
GGGTCCGCTGTTGGCAATACCGCCAGCTCTCTGCGCTGCACAGGGCTCCCCGCCACCCGGCCTGATAA
AGCTCGAAGACTGGGATACAAGGCTAAGCAAGGCTATGTCATTTACAGGATTCGTGTCCGCCGTGGTGGT
CGTAAACGCCAGTTCCTAAGGGTGCAACTTACGGCAAGCCTGTCCACCATGGTGTAAACCAGCTGAAAT
TTGCCGAAGCCTTCAGTCTGTTGCTGAGGAGAGAGCTGGGCGCCATTGTGGGGCTTTGAGAGTCCCTGAA
TTCTACTGGGTTGGTGAAGATTCCACATATAAATTTTTGAGGTTATCCTCATTGATCCGTTCCATAAA
GCTATCAGAAGAAATCCTGACACCCAGTGGATCACAAACCAGTCCACAAACACAGAGAAATGCGTGGGC
TGACATCTGTGGCCGAAAGAGCCGTGGCCTTGAAAGGGCCACAAGTTCACCACACTATTGGTGGTTC
CCGCCGTGCAGCCTGGAGAAGGCGCAATACTCTCCAGCTCCACCGTTACCGCT**TAA**

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_025586
Insert Size: 615 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_025586.3</u> , <u>NP_079862.1</u>
RefSeq Size:	1857 bp
RefSeq ORF:	615 bp
Locus ID:	66480
UniProt ID:	<u>Q9CZM2</u>
Cytogenetics:	14