

Product datasheet for MC210432

Gm17019 (NM_182957) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Gm17019 (NM_182957) Mouse Untagged Clone
Tag: Tag Free
Symbol: Gm17019
Synonyms: 4933431D05Rik; SPEER-4C
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >MC210432 representing NM_182957
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGGCCTTGTTTGCCCGACTCTGCAGACTTTTTCAGAGAGCCAATGTGGATGGGAGAGATACCAGAGAGG
GGAGGAAGGATGCTGACCTTCCATCTGAGAGGAATGAAGGACGAAGGAGGTGGACTTGCGAATGTGGAT
GGCTCTCAGACAGACATCATCCCTGTACCTGTCATAAGCAAGAAGCAGTTTGAGAAGGAAGAGAAAAAG
CTGATCAAAGAGATCCAGCTCACTACCGAGGAGACAAATGAGCTGAGAGATCGCCTGATCTACGTGACCG
AGGGATCCATGAACAAGAGGCCCTACCACAGGCAAAATCCCTTTATGAAAAATTGAAGTTAAAGGAGAA
GGAGATCATGACATTTCTGCACAACTTAGAGATGGAGAACATGGAGGCCAAGAGAACAAGCAGGAGCTC
AAGAAGGAGACACATTTCTATCGTAACTGCACAGCCGGCTCCTGATGGAGGAGAATCTTATAAAGAAGA
AGTCGATGACATTGCAGCAGGAGAGCAAGGAAGTACAGGCTGATTGGCCATCATTCACCAACGTTTGGT
GGAATTGAACCTGAGTGTAAAAGACGAACAGGAGAACAGCAATCTTGAGACTCCAGAATATCAGGTCTCA
GAAGCTGCAAGAGAGCTGGGATTGGCCACAGCTGAAGAAGACAGCATCCTGCAGAATGAGTTGCCAGGCC
AGGAGGCCCTGCTGAGCACCATCTCCAGCACCAATCTTCCCGGATGAATCTTCCATA**TAA**

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_182957
Insert Size: 768 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_182957.1</u> , <u>NP_892002.1</u>
RefSeq Size:	1171 bp
RefSeq ORF:	768 bp
Locus ID:	66773
Cytogenetics:	5 A2