

## Product datasheet for MC210446

### Serbp1 (NM\_025814) Mouse Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Serbp1 (NM\_025814) Mouse Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Serbp1  
**Synonyms:** 120009K13Rik; 9330147J08Rik; AL022786; Pairbp1  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**Fully Sequenced ORF:** >MC210446 representing NM\_025814  
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGATCGC**C

ATGCTGGGCACCTACAGGAAGGCTTCGGCTGCGTGGTACCAACCGATTCGACCAGCTATTTGACGACG  
 AATCGGACCTTTTCGAGGTACTGAAGGCAGCAGAGAACAAGAAAAAGAGCCGGCGGGGGCGGCTTGG  
 GGGCCCCGGGGCCAAGAGCGCGGCTCAGGCCGGCCAGACCACTCCAACCGCGGGGCAACAGTTG  
 CGTAAAGAGTCCCAGAAAGACCGCAAGAACCCGCTGCCCCAGCGTCGGCGTGGCCGACAAAAAGGAGG  
 AGACGCAGCCCGGTTGGCGCTTAAGAAAGAAGGAATAAGGCGAGTTGGAAGAAGACCCGATCAACAACT  
 ACAGGGTGATGGGAACTAATTGATAGGAGAGCAGAAAGGCGACCAACCCCGTGAAGAAGATTTGAAAAG  
 CCACTTGAAGAAAAAGGTGAAGGAGGTGAATTTTCAGTTGATAGACCGATTATCGAACGGCCTATCCGAG  
 GCCGAGGTGGTCTTGAAGGGGTCGAGGAGGCCGTGGACGTGGAATGGCCGAGGCGATGGATTTGATTC  
 TCGTGGCAAGCGTGAATTTGATAGGCATAGTGGAAGTGATAGATCTTCTTTTTCACATTACAGTGGCCTG  
 AAGCATGAGGACAAACCGGAGGTAGCGGCTCACAACCTGGGAACTGTCAAAGATGAATTAACAGAGT  
 CACCCAAATACATTAGAAACAAATATCTTATAATTGCAGTGATTTGGATCAATCAAATGTGACTGAGGA  
 AACACCTGAAGGTGAAGAGCACCTGTGGCAGATCTGAAAATAAGGAGAACGAAGTTGAAGAGTTAAG  
 GAAGAGGGTCCAAAAGAGATGACTCTGGATGAGTGGAAAGCTATTCAAAAATAAGACCGAGCAAAAGTAG  
 AATTTAATATCCGAAAACCAATGAAGGCGCCGATGGACAATGGAAAAAGGGATTTGTTCTGCATAAATC  
 AAAAAGTGAAGAGGCTCATGCTGAAGATTCAGTTATGGACCATCATTTCCGGAAGCCAGCAAATGATATA  
 ACATCTCAACTGGAGATCAATTTTGGAGACTTAGGCCGCCAGGACGTGGTGGCAGAGGAGGACGTGGT  
 GCGTGGCGTGGTGGACGTCAAACCGTGGCAGCAGGACTGATAAGTCAAGTGCTTCTGCTCCTGATGT  
 AGATGACCCAGAGGCCTCCAGCTCTGGCT**TAA**

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
 ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

<b>Restriction Sites:</b>	Sgfl-Mlul
<b>ACCN:</b>	NM_025814
<b>Insert Size:</b>	1224 bp
<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_025814.2</a> , <a href="#">NP_080090.2</a>
<b>RefSeq Size:</b>	6685 bp
<b>RefSeq ORF:</b>	1224 bp
<b>Locus ID:</b>	66870
<b>UniProt ID:</b>	<a href="#">Q9CY58</a>
<b>Cytogenetics:</b>	6 C1
<b>Gene Summary:</b>	May play a role in the regulation of mRNA stability. Binds to the 3'-most 134 nt of the SERPINE1/PAI1 mRNA, a region which confers cyclic nucleotide regulation of message decay. Seems to play a role in PML-nuclear bodies formation.[UniProtKB/Swiss-Prot Function] Transcript Variant: This variant (1) encodes the longest isoform (1).