

Product datasheet for **MC210462**

Trappc13 (NM_001093759) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Trappc13 (NM_001093759) Mouse Untagged Clone
Tag: Tag Free
Symbol: Trappc13
Synonyms: 2410002O22Rik; 2610524F24Rik; AI451892; AU017229
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >MC210462 representing NM_001093759
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**GCGATCGCC**

ATGGAAGTGAATCCCCCAAACAGGAGCACCTGCTGGCGCTAAAAGTGATGCGATTGACTAAACCTACTT
 TGTTACCAATATTCCTGTAACTGTGAAGAGAAGGACTTACCTGGTGATCTCTTAAACCAAGTAAATGAA
 AGATGATCCTTCAACTGTAATGGTGCAGAGATTTAATGTTGGGAGAAATGTTGACTTTACCCAAAAT
 TTCGGGAATATATTTTTGGGAGAGACCTTTCCAGTTATATCAGTGTTCAATGACAGCAATCAAGTTG
 TAAAAGATATATTGGTGAAGCGGATCTGCAGACGAGTCTCAGCGTTTAAACCTTTACGCTCCAATGC
 TGCTGTGGCTGAACTTAAACCCGACTGCTGATTGATGACGTCATACCCATGAAGTAAAGGAAATTGGG
 ACTCACATCTTGGTGTGTGCTGTGAGTTACACAACGCAGGGTGGAGAAAAAATGACTTTAGGAAATTC
 TCAAATTCAGGTTCTCAAACATTGGATGTGAAAACCAATTTACAATGCAGAGAGTGACCTCAGTTC
 TGTGACTGATGAAGTATTTCTGAAGCACAAATTCAGAATATTACAACCTCACCAATGTTTATGGAGAAA
 GTTTCACCTAGAGCCATCAATAATGTACAATGAACAGAATTAACCTCGGTGACCCAGGCTGGAGAATGCA
 TATCTACATTTGGGTCCAGAGGATATTTGCAGCCAATGGATACACGCCAATACTTATACTGCCTAAAGCC
 CAAGAAGGAATTTGCAGAAAAGCCGGCATAATTAAGGGAGTGACAGTAATTTGGGAACTGGATATAGTG
 TGGAAAACAAATCTAGGTGAAAGGGGACGGTTACAGACCAGCCAGCTTCAAAGGATGGCTCCAGGTTATG
 GGGATGTTAGGCTGTCTTTAGAGGCCATCCCAGACACTGTAACCTAGAAGAACCTTTCCATATTACCTG
 TAAATCACAACTGCAGTGAGAGGATGATGGACCTCGTGCTGGAGATGTGCAACACCAATTCATCCAC
 TGGTGTGGGATCTCAGGACGACAGCTGGGAACTGCACCCAGCTCCTCGCTCTGCCTGGCCCTGACGC
 TCCTGTCTCAGTACAGGGGCTCAAAGTGTCTCTGGCTTAAGACTCACTGACACATTTCTAAAGAGAAC
 ATACGAATATGATGACATCGCGCAAGTCTCGTAGTCTCAGCAGCCGTGGAAGTGAAGCC**TAA**

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA



[View online >](#)

Restriction Sites:	Sgfl-Mlul
ACCN:	NM_001093759
Insert Size:	1254 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001093759.1</u> , <u>NP_001087228.1</u>
RefSeq Size:	3560 bp
RefSeq ORF:	1254 bp
Locus ID:	66975
UniProt ID:	<u>Q3TIR1</u>
Cytogenetics:	13 D1