

## Product datasheet for MC210549

### Nts (NM\_024435) Mouse Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Nts (NM_024435) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Nts
Synonyms:	5033428E16Rik; NMN-125; NN; NT; NT/N; NTS1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>MC210549 representing NM_024435 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGAGAGGAATGAATCTCCAGCTGGTGTGCCTGACTCTCCTGGCTTTCAGCTCCTGGAGTCTGTGCTCAG  
ATTCAGAAGAAGATGTGAGAGCCCTGGAGGCAGATCTATTGACAAACATGCATACATCCAAGATCAGCAA  
AGCAAGTCCTCCGCTTTGGAAAATGACCTTGCTAAATGTTTGACGCTCATAAATAACGTGAACAGCCCG  
GCCGAGGAAGCAGGAGACATGCATGATGACGACCTTGTGGCAAAGGAACTTCCCTTGTCTGGATG  
GTTTTAGCTTGAAGCAATGCTGACCATCTCCAGCTCCAGAAAATCTGCCGAGCAGGGCCTTTCAACA  
CTGGGAGATAATCCAGGAAGATATCCTTGATAACGTCAATGATAAAAACGAGAAGGAAGAAGTGATAAAG  
AGAAAAATCCCTTATATTCTGAAAAGGCAGCTGTATGAAAATAAACCCAGAAGGCCCTACATTCTCAAGA  
GGGTTTCTACTACTACTGA

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	Sgfl-MluI
ACCN:	NM_024435
Insert Size:	510 bp



[View online »](#)

**OTI Disclaimer:** Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at [custsupport@origene.com](mailto:custsupport@origene.com) or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_024435.2](#), [NP\\_077755.1](#)

**RefSeq Size:** 1237 bp

**RefSeq ORF:** 510 bp

**Locus ID:** 67405

**UniProt ID:** [Q9D3P9](#)

**Cytogenetics:** 10 D1

**Gene Summary:** Neurotensin may play an endocrine or paracrine role in the regulation of fat metabolism. It causes contraction of smooth muscle (By similarity).[UniProtKB/Swiss-Prot Function]