

## Product datasheet for **MC210596**

### Trmt112 (NM\_001166370) Mouse Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Trmt112 (NM_001166370) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Trmt112
Synonyms:	0610038D11Rik; Trm112p
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>MC210596 representing NM_001166370 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGCATCGCC**

ATGAAACTGCTCACCCACAATCTCCTAAGCTCCCATGTGCGGGGCGTGGGCACGCGTGGCTTCCCGCTGC  
GCCTGCAGGCCACTGAGGTCCGCATCAACCCCGTGGAGTTCAACCCGGAGTTCGTGGCGGGATGATCCC  
CAAAGTGGAGTGGCGGCGCTTGTGCAGGCAGCGGACACCTTAACTTGGCCGAGGTACCCAAAGAGCCG  
ACTGAAGGCTATGAGCACGACGAGACATTTTTGAGGAAGATGCACCACGTGTTGCTTGAAGTTGATGATC  
TGGAGGGCACCTGCAATGCCCAAGATCTGGTCGCTTTTCCCATCAGCCGCGGGATCCCAATATGCT  
GTTGAATGATGAGGAGACTGAGACG**TAA**

AG**CGGACCG**ACGCGTACGCGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCC  
TGGATTACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	SgfI-RsrII
ACCN:	NM_001166370
Insert Size:	378 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_001166370.1</a> , <a href="#">NP_001159842.1</a>
<b>RefSeq Size:</b>	1047 bp
<b>RefSeq ORF:</b>	378 bp
<b>Locus ID:</b>	67674
<b>UniProt ID:</b>	<a href="#">Q9DCG9</a>
<b>Cytogenetics:</b>	19 A
<b>Gene Summary:</b>	<p>Acts as an activator of both rRNA/tRNA and protein methyltransferases (PubMed:26797129). Together with methyltransferase BUD23, methylates the N(7) position of a guanine in 18S rRNA (By similarity). The heterodimer with HEMK2/N6AMT1 catalyzes N5-methylation of ETF1 on 'Gln-185', using S-adenosyl L-methionine as methyl donor (PubMed:20606008, PubMed:26797129). The heterodimer with ALKBH8 catalyzes the methylation of 5-carboxymethyl uridine to 5-methylcarboxymethyl uridine at the wobble position of the anticodon loop in target tRNA species (By similarity). Involved in the pre-rRNA processing steps leading to small-subunit rRNA production (By similarity).[UniProtKB/Swiss-Prot Function]</p> <p>Transcript Variant: This variant (1) represents the longer and more frequently occurring transcript. Both variants encode the same protein.</p>