

Product datasheet for **MC210627**

Snrnp48 (NM_026382) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Snrnp48 (NM_026382) Mouse Untagged Clone
Tag: Tag Free
Symbol: Snrnp48
Synonyms: 1110050F08Rik; 6530403A03Rik; C6orf151
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >MC210627 representing NM_026382
Red=Cloning site **Blue**=ORF **Orange**=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGGAGACCGAACCTCCTCCTCTGGAAGAGCGGCGTCGGCTCCAGGAGGAGCTGAGCGAGTTTGTGGAGA
GCTGCTGCCGGACGCTGGAGGAGGTGACGGCGTCCCTGGGATGGAGCCTGGATCAGCTGGACCCCGGGA
CGAGGCGGAGGCCGAGGATGAAATCGCCATATGTCCTTATGATTCCAATCATCGCATGCCAAATCATCA
CTAACAAAGCACATGGAATCTTGACAGATTGAGGAACTGGGATATACCAAGAGGAAGAGAATGAAATGT
ACAATCCTACTTTTTTCTATGAAAATTTGAAGATTCTTCAGTTACTTTGAATAAGGATTCACAGTTCCA
GATAATTAACAAGCTAGGACTACAGCTGGGAAAGATGGTATTGTTATAGCCAGAGAATGTATTCTTCA
GTGCTGTTGAAGTTCTCTGAATCACAAACGATCTGTTTGTGACCTCACACAAACGGATCGCCTTGCC
TGTATGATTTTGTGATTGAGGAGACAAAGAAAAGCGCTCTGGTCTCAAGTCATTGAAAATGACAGTGA
CCTGTTTGTAGACTTGGCTGCTAAAGTCAATCAAGATAATAGCCGAAAAGTCCAAAATCTTACCTGGAA
ATCCTGGCAGAAGTGAGAGATTATAAAGAAGACGCCAGTCTACAGAGCTAAGAATGTTACATAACCA
AGAAATCATACTACTGAGGTGATCCGGGATGTGATCAAGGTGCACATGGAGGAGCTGAGCAGCCACTGGCA
GGAAGAGCAGGGGAGAGCAGGAGACGCTGCTGAGAAGAATGAAGAAAGGCGGTACGCTTACGTAGATTCA
CGGCAGTCTGGTGAAGTTATTTGGATGTTGAAAGTTCACGCATAGAAGGGCTCGGAGTAGGAGCCAC
ACAAACGAAAAGAAATAAAGATAAAGTTCTGAGTCAAGAAGAAGAAAAGAGAGGGATGGAGAAAAGACA
TCATAGTCATAAAGAAGAAAGCAAAAATATA**A**

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_026382



[View online »](#)

Insert Size:	1014 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_026382.2 , NP_080658.2
RefSeq Size:	1789 bp
RefSeq ORF:	1014 bp
Locus ID:	67797
UniProt ID:	Q9D361
Cytogenetics:	13 A3.3
Gene Summary:	Likely involved in U12-type 5' splice site recognition.[UniProtKB/Swiss-Prot Function]