

Product datasheet for **MC210766**

1110012L19Rik (NM_026787) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: 1110012L19Rik (NM_026787) Mouse Untagged Clone
Tag: Tag Free
Symbol: 1110012L19Rik
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >MC210766 representing NM_026787
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGAAGTTTCCCTGCCTATCCTTCCGGCAGCCCTACGCTGGTTAATCTTAAATGGAGTCAAGACTCTTG
AGACACGTTGGCGTCCCCTGTTGAGCAGCGTTCAGAAATACACCATTGCCATCCACATTGCTCACAAGGA
CTGGGAAGATGATGAATGGCAGGAAGTGCATGGAGAGGTTGGGATGACCTGGACACAGATTCAGACC
TTACTACAGGCAGGAGAAAAGTACGGCAGAGCGTGATAGCTGGACTTATTGACATTGGAGAGACTTTTC
AGTGTCCAGAGACCTTAACTGCTGAGGAGGCTGTGAACTGGAACTCAAGCTGTGCTAACCAATCTGCA
GCTGAAGTACCTGACTCAAGTTTCAAACCCAGGTGGTTACTGGAGCCATACCTAGGAAAGGGGAAAG
GATATATTCCAGGTAGACATCCCAGAGCACTTGATCCCACTGGAGAAAGAA**TGA**

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_026787
Insert Size: 474 bp
OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).



Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_026787.2](#), [NP_081063.1](#)

RefSeq Size: 1057 bp

RefSeq ORF: 474 bp

Locus ID: 68618

UniProt ID: [Q9D1F3](#)

Cytogenetics: X A7.2

Gene Summary: May play a role in cell protection during the inflammatory response. In epithelial cells, negatively regulates IL6 production and apoptosis through the regulation of MT2A expression. [UniProtKB/Swiss-Prot Function]