

Product datasheet for **MC210770**

Scgb3a1 (NM_054037) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Scgb3a1 (NM_054037) Mouse Untagged Clone
Tag: Tag Free
Symbol: Scgb3a1
Synonyms: HIN; HIN-1; LuLe; LuLeu2; Pnsp2; UGR; UGRP2; Utgrp2
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Fully Sequenced ORF: >MC210770 representing NM_054037
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGAAGCTTACCACCACCTTTCTAGTGCTCTGTGTGGCTCTGCTCAGTGACTCTGGTGTGCTTTCTTCA
TGGACTCATTGGCCAAGCCTGCGGTAGAACCCGTGGCCGCCCTTGCTCCAGCTGCAGAGGCTGTGGCAGG
GGCTGTGCCTAGCCTACCATTAGCCACTTGGCCATCCTGAGGTTTCATCCTGGCCAGCATGGGCATCCCA
TTGGATCCTCTCATAGAGGGATCCAGGAAGTGTGTACCGAGCTGGGCCCTGAGGCTGTAGGAGCTGTGA
AGTCACTGCTGGGGTCTGACAATGTTTCGGT**TGA**

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Chromatograms: https://cdn.origene.com/chromatograms/ja2795_a03.zip
Restriction Sites: SgfI-MluI
ACCN: NM_054037
Insert Size: 315 bp



[View online »](#)

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_054037.2](#), [NP_473378.1](#)

RefSeq Size: 540 bp

RefSeq ORF: 315 bp

Locus ID: 68662

UniProt ID: [Q920D7](#)

Cytogenetics: 11 29.75 cM

Gene Summary: This gene encodes a small secreted protein that may function in the inhibition of cell growth. Expression of the related gene in humans is downregulated via DNA methylation in breast cancers and other tumor types. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene. [provided by RefSeq, Jan 2013]
 Transcript Variant: This variant (1, also known as type A) represents the shorter transcript and encodes the longer isoform (a).