

Product datasheet for **MC210830**

Tamm41 (NM_026894) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Tamm41 (NM_026894) Mouse Untagged Clone
Tag: Tag Free
Symbol: Tamm41
Synonyms: 1500001M20Rik; TAM41
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >MC210830 representing NM_026894
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGGCGCTGCAGGCGCTGCACAGCTCGGGCGTGGGCTGCGCAGGATCCTGGCTCACTTCCCGGAGGACC
TGAGTTTGGCCTTCGCCTACGGCTCCGCGGTGTATCGCCAGGCAGGACCCAGTCCCACCAGGAGAATCC
CATGTTAGACTTAGTGTTCACAGTAGACGACCCTGTTGCCTGGCAGCGATGAACCTGAAGAAGAATTGG
AGCCACTACTCCTTCCCTCAAACCTCTGGGGCCAGGATCATCTCCTCCATCCAGAATAATTACGGTGCCG
GAGTCTACTTCAACCCCTGATCAGGTGTGACGGGAAGCTTATCAAGTACGGGGTTATTAGCACTGGCAC
ACTGATTGAGGACCTCCTGAACCTGGAATAACCTATATATCGCCGGACGCCTCCAAAAGCCGGTGAAGATC
GTATCGATGAATGAGAACATGGCTCTTAGAGCCGCCCTGGATAAAGAACCTGAGGAGTGCCGTGACCACCG
CCTGCCTCATGCTACCTGAAAGCTTCTCCGAGGAAGACCTGTTTCATAGAGATAGCTGGGCTCTCCTACTC
AGGTGACTTCCGCATGGTATTGGCGAGGAAAAATCAAAGTGTGAACATCGTGAAGCCCAACGTAGGC
CATTTCCGAGAGCTCTATGAGAGTATATTACAGAAGGACCCCAAGTGGTGTATAAGATGCATCAAGGCC
AACTGGAGATAGATAAAAGCCAGAAGGACAGTTCACCTCAGCTAATGACACTGCCAGAACCCTACAGCA
GCAAATCAATCACATCATGGACCCTCTGGGAGAAACAGAGACGTGGAAGAGACTCTGCTGCAGGTGGCT
CAGGACCCCGACTGTGGAGATGTGGTGCAGTACTCATCAATCGTCAGACCTTCCAGCATCAGGC
AGAGCACTAAAGGCCTTTTACTGCTGGTATGAAGAAGTCTGTGATTTATAGTTCACGGAAACTGAACAA
AATGTGGAAGGGTGGATGAGAAAAGCGTCC**TGA**

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_026894



[View online »](#)

Insert Size:	1014 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_026894.1</u> , <u>NP_081170.1</u>
RefSeq Size:	1183 bp
RefSeq ORF:	1014 bp
Locus ID:	68971
Cytogenetics:	6 E3