

Product datasheet for MC210869

Mcemp1 (NM_026985) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Mcemp1 (NM_026985) Mouse Untagged Clone
Tag: Tag Free
Symbol: Mcemp1
Synonyms: 1810033B17Rik
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Fully Sequenced ORF: >MC210869 representing NM_026985
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTGTAAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGCATGCATCAGCCTCCCAGGATAAGAACCGGAGGAAGCCAGGTCATGATGAAGGTGCTCACAATCCTG
 ACTACGAGAATATAACCTTGGCCTTCAGAAACAAGGACCAACTCAAACCTCAGCCAATCAACACCCACAAA
 ACAAGCCAAGTTCAAGACATCCCTGGACCCAGCTGAGTCCCCGCCTTGGTTGTACAGAACCATTATGATG
 TTGTATGTTCTCCTTGTCTCGTCTTTTATCCTGCATCGTCTCTCTGCTTTGGTCTTGGTGAAAAATT
 CTGAGATGTCCAAGGAGCTGTGGACCTTGAAGCAGAGCTTTCGAATGTTTCAGACACGGTGTGGAATAT
 CCGGGAGCTCCAGAATCAGCAAACGAGGATTTGGGAAGCTGCCAGGGGGACATCAAGGAGGTCAAGAAG
 ACCCTTGGCACAGTCATGAGTAGCATCCAGACTGGAAACGACCGCTGAAGACTGTGCCGCGCAGATATAA
 CCCAAATCAAGAAAACCTTTGAGGCGCTAGAAAAGAAGGCACAGCCTCAGCCCAGTACA**TAA**

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_026985
Insert Size: 552 bp
OTI Disclaimer:

Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).


[View online »](#)

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	NM_026985.1 , NP_081261.1
RefSeq Size:	1132 bp
RefSeq ORF:	552 bp
Locus ID:	69189
UniProt ID:	Q9D8U6
Cytogenetics:	8 1.92 cM