

Product datasheet for MC210893

Izumo3 (NM_027034) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Izumo3 (NM_027034) Mouse Untagged Clone
Tag: Tag Free
Symbol: Izumo3
Synonyms: 1700011H22Rik
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Fully Sequenced ORF: >MC210893 representing NM_027034

Red=Cloning site Blue=ORF Orange=Stop codon

TTTGTAAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGGAGACCTATGGGTACTTCTGTTTTTCATCTCTATCCCTAGCAGCCTTTCATGGGGTCAGAGGCTGTT
TGGAGTGTGACCCCAAATTTACAGAGGATGTTAGGACCTTGTGGGAAACCTGGTACCTCTAGAAGTTCC
TGATAGAAATCACTACTGGACAGGCAACATAAGGAGATAACCCATATAAGCTCCAAGGTCTCCACAAG
GATAAGATGCTTCGGTTGTTGGCTGTTGAAATGTGATCAAATTGAGAGAATGGCTGAAGGATGAATTTT
ATAGGCTGGGCAACGAAACATGAAAGGAGTCTTTATCATTCAAGGCAATCTTCTGGAGATCCGCCAAAA
CTTGGAGAAAAAAGTACACAAATAATAAAGAACTTCTCTGAAGTTGCTTGTTCAGAAGATTGCATTGTA
ATTGAAGTCTGTCTCGATTGTTGAACTGTCTTCGGATGACCACTCGGTGTTTCAAAGGAGACTATT
GTCAAGATGAGGATCCAAAAACCGCTGAAATCGAGAGATTTCCCTCTATCTGATATTTATAGCAGAAGC
TGTCATTTTGGCGAGTGCTGTGCTACTATTCCATGTTTGCATCTCACATCGGAGGAAAAATGAAGGCAATA
CGAAGAACTTTGAAGACATATTTGGAGAAGAACTAGAAGAACTAGTGGAGATGATATACAAAGATGATG
GTGAAAAAATTTGGAGTTTGGGAGATCAAATTCAACTCTCTGACAGGAGAGCCTACATGTGCTGAGTC
AGAAATGCAAACTGGGACT**TGA**

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_027034
Insert Size: 792 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_027034.1, NP_081310.1</u>
RefSeq Size:	1263 bp
RefSeq ORF:	792 bp
Locus ID:	69314
UniProt ID:	<u>A6PWW3</u>
Cytogenetics:	4 C5