

Product datasheet for **MC210896**

1700012B07Rik (NM_001162428) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	1700012B07Rik (NM_001162428) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	1700012B07Rik
Synonyms:	1700069D21Rik
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin
Fully Sequenced ORF:	>MC210896 representing NM_001162428 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGTCCCATTTCATAATTGGAAGAGCGCCAAACGAGCCAGAGCTTGACGTACAGGACAGCATCAGGG
TTTCTCATGAGAGAAAAGTGAGAACAAAATCATGAAGAGCTTCCAGAGAGCCACCACTCATCCCAGG
TCAGGACAAGGAGGGAGTCTGGCTTCTGGGAGCAGAAGTGAAGAGGCAGAAGAAGACCCGGCTGCAGAG
GAGGGACCTTTTCTCAGTAAAACAGAAGTACAAAAGATACGCAAGGAGCATGTGACAGCAAAGTCTGTTG
GACAAGGTCTGTAATTGTCCCCTCCAACAACAAAGAGCATCCGAATATACGACGCGAGCCCTTTTCTTC
TTACCCGTTGACAGGCTGAGAAGATGAGTCTCACAGAAGGGCACGTATGGCCCGTCTGCCACGACAGGGC
CCTTTCCGGCACCTTATGGACATGAAGGCAGAGAAGAGGCTGGCGGGGTGGAAGGAACGCCGTAGCCGCT
TTGGTGACATCAACGTGCACAAATGCTATCAGTTTGGGCAGGTCTTCTCTCTTTCCAGGCTCTGGCCAC
CACAGAGATGCCAAGGCTGGTGTGCCAGAGACTTGCCCATGAGTTCTCACAGGCGGAGAATAAACCTT
TCATCCCTCGGTGACATGAACCTCCAGGATTCACCCCTGTCTCCAGAGCTGACATTGGGAAACGATGATG
GCCATCATGAGAAGGACAAGCCAGTCCATCATTCAATCATCAGGGTCTCTCCACGTGATGGACGTTAC
AGAAGCCTCAAACCTCTGCTTGCATCTCATCACTGCGCCTGGAAGACTATTCCATCTCTGA

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	Sgfl-MluI
ACCN:	NM_001162428
Insert Size:	831 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001162428.1</u> , <u>NP_001155900.1</u>
RefSeq Size:	1395 bp
RefSeq ORF:	831 bp
Locus ID:	69324
Cytogenetics:	11 E1