

## Product datasheet for **MC210917**

### Zpbp2 (NM\_027061) Mouse Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Zpbp2 (NM\_027061) Mouse Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Zpbp2  
**Synonyms:** 1700017D11Rik; 2610022C02Rik; ZPBPL  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**Fully Sequenced ORF:** >MC210917 representing NM\_027061  
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGGATCGCC**

ATGCTTGGCTGGGCCCTACTCTACGCGGTGCTCTGGTCTCTCGCAGGAGTTGGATCGCAACGTTCTTCTCT  
TATTCAATATCAAAGCCTTTGTTACGGCACGGTAGGGCACCCAGTAAAAATATATGTAAGTTACATCA  
GGCTAGTCTGTCTTATTTGTATGGATATTGATCGTGCCAATAAAGAAACAGTGGACCCACCTATTTA  
TGGACTGGGCCAAATGAAAACACATTAAGGGAATAGCCAGATAAATATAACAAACACAGGAGAACTGG  
TACTGAAGGATTTCTGGAGCCGTTATCTGGACTTTACTCATGTATGCTTTCCTATAAGACTATCAAAGC  
AGAAACCAAGAAGAACTATGATAAAGAAGAAATATGACTTTCTGGTCTTTGCATATCGGGAACCTGAT  
TATTCTTACCACATGGCCGTGCGTTTTACCACAGGATCTTGTGTAGGAAGACACAATGACCTGCTCTTTA  
GAGTGTGAAGAAAACTTGGATAACTTAATCTCTGATTTATTGTGCCATGTCATAGGCCATCATAAA  
ATGCCATTTTGTAACTTCCAGAGAGGGACTTCTGTATGAGCTGTTTCATAGCCTTCAAGTGAATCCT  
TTCGCACCGGGCTGGAGAAGTTTGTGCAACCGCTCTGCTGACTGTGAAGACATCACGAATCACACGTCC  
TCAAGGCGAGAGATCGGATGGAAGAATTTTTCGGAAACAAGCACACATTCTACACCTCACTTCAACAG  
AACGGTCCCAGCTATGCACTTTGTGGACCACAGTTTTCAAGTAACACGGATAGACAATTGTCGACCAGGC  
TTTGGGAAAAACGAAGTCTGCACAGCAATTGTGCTACCTGTTGCGTGGTCTGTAGTCTGGAACGTTCA  
GTCCCATGTGGATGTCACCTGTCAGATCTGTGTTCTGTCCATATCTATGGAGCTAAAGCTTGCCATA  
A

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** SgfI-MluI  
**ACCN:** NM\_027061



[View online »](#)

|                               |   |
|-------------------------------|---|
| <b>Insert Size:</b>           | 981 bp  |
| <b>OTI Disclaimer:</b>        | Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).  |
| <b>Components:</b>            | The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).  |
| <b>Reconstitution Method:</b> | <ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol> |
| <b>RefSeq:</b>                | <a href="#">NM_027061.2</a> , <a href="#">NP_081337.2</a>   |
| <b>RefSeq Size:</b>           | 1519 bp   |
| <b>RefSeq ORF:</b>            | 981 bp  |
| <b>Locus ID:</b>              | 69376   |
| <b>UniProt ID:</b>            | <a href="#">Q6X786</a>  |
| <b>Cytogenetics:</b>          | 11 D  |
| <b>Gene Summary:</b>          | Is implicated in sperm-oocyte interaction during fertilization.[UniProtKB/Swiss-Prot Function]<br>Transcript Variant: This variant (1) encodes the longest isoform (1).   |