

Product datasheet for **MC210964**

Vkorc111 (NM_027121) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Vkorc111 (NM_027121) Mouse Untagged Clone
Tag: Tag Free
Symbol: Vkorc111
Synonyms: 2310024K08Rik
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >NM_027121.3
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGGCGGCGCCCGTCTGCTGAGAGTGTCCGGTCCCGCTTGGGAACGGGTGGCCCGGTATGCAGTGTGCC
 CCGCCGGGATCCTGCTCTCCATCTACGCCTACCACGTGGAGCGGGAGAAGGAGAGGGACCCGGAGCACCG
 GGCCCTCTGCGACCTGGGGCCCTGGGTGAAGTGTCTCCCGCCCTGGCCCTCAGATGGGGTGCAGGATTT
 GGTCTTTGGGTTCCATTTTGGAAAAGATGGTGTATTAACAGCCAAACAGTGTCTTTGGACTTATAT
 TTTATATACTACAGTTATTACTTGGCATGACAGCCAGCGAGTTGCAGCTCTGGTCTCATGACCTCCTC
 CATTGTGTCTGTGGTGGCTCTTTGTACCTGGCTACATTCTGTACTTTGTGCTGAAAGAGTTTTGCATC
 ATCTGCGTCACCACATATGTGCTGAACCTCCTCCTCATCAATTACAAACGACTAGTTTATTTGA
 ATGAGGCTGGAAGCGACAGCTGCAGCCTAAGGAAGACTGA

ACGCGTACGCGGCGCTCGAGCAGAACTCATCTCTGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: Sgfl-MluI
ACCN: NM_027121
Insert Size: 531 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_027121.3, NP_081397.1</u>
RefSeq Size:	1113 bp
RefSeq ORF:	531 bp
Locus ID:	69568
UniProt ID:	<u>Q6TEK5</u>
Cytogenetics:	5 G1.3
Gene Summary:	<p>Involved in vitamin K metabolism. Can reduce inactive vitamin K 2,3-epoxide to active vitamin K (in vitro), and may contribute to vitamin K-mediated protection against oxidative stress. Plays a role in vitamin K-dependent gamma-carboxylation of Glu residues in target proteins. [UniProtKB/Swiss-Prot Function]</p> <p>Transcript Variant: This variant (1) represents the longest transcript and encodes the longest isoform (1).</p>