

## Product datasheet for **MC211086**

### Zdhhc24 (NM\_001168517) Mouse Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Zdhhc24 (NM\_001168517) Mouse Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Zdhhc24  
**Synonyms:** 5730496N17Rik; Leng4  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**Fully Sequenced ORF:** >MC211086 representing NM\_001168517  
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**GCGATCGCC**

**ATGCTACTCACAGCAAAGTGTCTCTGGCCAGTTGCCCTGGCCTTTGGTGGACACGTGTGTGGCAG**  
**GTGCACGTGTGTGGTCTGGGCTGCTTTCCATGGGATGCTGCTGCTTCGGGGACAGACCACCTGGGA**  
**GTGGGCTCGGGTCACTACTGCTATGACCTGGCACCTGCCACAACCTGCAGGCAGCCCTGGGGCCTCGC**  
**TGGGCCCTAGTCTGGTTCTGGCCTTTCTAGCCTCTCCTTGGCTGGAGATGGGATCTCTTCCAGACTC**  
**CAGGTGATGTGGTCTTGTGACTTC**TGA****

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** SgfI-MluI

**ACCN:** NM\_001168517

**Insert Size:** 309 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).



[View online »](#)

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** NM\_001168517.1, NP\_001161989.1

**RefSeq Size:** 2316 bp

**RefSeq ORF:** 309 bp

**Locus ID:** 70605

**Cytogenetics:** 19 A

**Gene Summary:** Probable palmitoyltransferase that could catalyze the addition of palmitate onto various protein substrates.[UniProtKB/Swiss-Prot Function]  
Transcript Variant: This variant (3) differs in the 5' UTR, lacks a portion of the 5' coding region, and uses a downstream start codon, compared to variant 1. The resulting isoform (b) is shorter at the N-terminus, compared to isoform a. Both variants 2 and 3 encode the same isoform. Sequence Note: This RefSeq record was created from genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.