

## Product datasheet for **MC211090**

### Oip5 (NM\_001042653) Mouse Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Oip5 (NM\_001042653) Mouse Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Oip5  
**Synonyms:** 5730547N13Rik; Lint-25; Ms18b  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**Fully Sequenced ORF:** >MC211090 representing NM\_001042653  
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**GCGATCGCC**

ATGGCGACTCTCTCGCGCCGAGTCTTCGTTCTTGGAGACCCAGTCACGTCGTGACTCTTGTGATACAC  
CTTTTTCTTACCTCCTCCATGGAGTGGGACACTCAGGTGGTGAAGGGTTCTCCCTCTCGGGTCTCC  
CGAGTCTAAGGTGAAGGGCTCACCTGTCGACCTGCGGCTGCCTGCGTGGTGGAGCCGGAGAGATGTGCC  
GTGTTTCATTGCGCGCGGTGCTACGCGGTGCTCGGCGACACCCTGCACCTCGCCTGGGACCTGTCGCGGT  
CCCTCGGCGCCCTGGCCTTCTCCAAAGTCACAAACAACGTCGTTTTGTTAGAGCCCTTTTTGGTTGGCAT  
AGAAGTTTTCTCAAGTCCAGCACATACAACCTTTTATTCTGTAATTCCTGTGGGACTCCTGTTGGTTTC  
CATTTATATTCTACCCATGCTGCTATGGCTGCCTTGAGAGGCCATTTCTGCCTTCTAGTGATAAAATGC  
TGTGCTATCTGTTAAAAACAAATGCCATAGTAAATACATCTGAAATGGATTTTCATAATGTACCTTTACC  
AGAAAAGATTGCAGAGCTGAAAGAGAAGATAATGCTAATGCATACTCGTCTAAACTACTGACGGGGCTT  
CTGAAGGGAAAGTCTCCTCATCAGTTTAAAGCAAGAAAACCAACAAGCAAGGAAACAACACATCTTGGGGT  
TAACAGCATCACAAAGATCCTG**TAA**

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** Sgfl-MluI  
**ACCN:** NM\_001042653  
**Insert Size:** 726 bp



[View online »](#)

<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">NM_001042653.1</a></u> , <u><a href="#">NP_001036118.1</a></u>
<b>RefSeq Size:</b>	1557 bp
<b>RefSeq ORF:</b>	726 bp
<b>Locus ID:</b>	70645
<b>UniProt ID:</b>	<u><a href="#">A2AQ14</a></u>
<b>Cytogenetics:</b>	2 E5
<b>Gene Summary:</b>	Required for recruitment of CENPA to centromeres and normal chromosome segregation during mitosis.[UniProtKB/Swiss-Prot Function]