

Product datasheet for MC211117

4921539E11Rik (NM_001163494) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: 4921539E11Rik (NM_001163494) Mouse Untagged Clone
Tag: Tag Free
Symbol: 4921539E11Rik
Synonyms: RP23-47A6.2
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >MC211117 representing NM_001163494
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGATAAACTATGAAGAAAGTAAATTAATGGACCACAAAGAAGAAAATCTAAAATCAAGACCTATATCTC
TTCGTTATCTTAAGGATAAAGATGAGACAGAATACGCCAATCCCCTGCCTTTCCACAGCTTTGGTCTAA
ACACCTGTGTA AAAAATCAGCGGAGTCTTCAGTCCGACTCCCAAACCTACCAATGAATCCAATGCTAGT
AAGAAGGAAAACGTCTCTCCTCATTACAGATCAACATGAAAGCAGAATAAGAAGTCCATGCATTCTA
CAGACCACTCAGCAGATTCTAGTACAAGCAGGGGAAATGCCCTCCAAAGGGTATTACAAAAGAAAGTGA
GTTAACCGAAATGATGAAGCCAGAAAGCCCATCCAGTCAAACAGAGCATCATGCTTCCCTGGACTGT
GAGGACCTACTGAAAAATCCAAAATAAAGACAATTGACCTTCGTCCAGCAGTGACAGTACATACTTCCA
TGAAACAATCCCACGCAATCCAATAATTTCTATGACTGAATATGTACAAATGCTGTTCCCTGACCAA
AAGTTTACCCCTACGCTATGAAATGTACAGAAAGGAATATTGTATTAGAAAAAATTATGAAGTACTC
AAAGTTTGTAGTATGAACCAAGTGTCTTCTGAACCTATCCAGCAAAAACATTTGCAAGTCTTCT
CTGCTGAATATGCCAAAAGAGTATAAATGAGAAACGAAAGAAGAAACATGACCGTCTGGTTTCAAAGAA
AATATCTCCGAACACTCGTTATAATTTATCCAAACATTTCCAGCTTAAGCAAAAATTTGTTGGTTAC
TTTGATAAAGATGTTACTCAAGGAAAGAGTTAATAGGCAAACAGATTTGAAAGGTTTTCCAAAACAAGC
CACCACCACACGCAAACCTACTACCTTACCTATCAAGTATGACTCCAAGCCTCTGAAAAATATATTTGA
AATACATAAACTAAACAACATGACCCCATTTGGATAACCTGCTAGGCTTGAGAGCA**TAA**

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_001163494



[View online »](#)

Insert Size:	1038 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001163494.1</u> , <u>NP_001156966.1</u>
RefSeq Size:	2356 bp
RefSeq ORF:	1038 bp
Locus ID:	70941
UniProt ID:	<u>Q9D5Q8</u>
Cytogenetics:	4 C6