

Product datasheet for MC211148

Golga7b (NM_001141983) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids

Tag: Tag Free

Symbol: Golga7b

Synonyms: 4933417O08Rik; A1839934

Mammalian Cell Selection: Neomycin

Vector: pCMV6-Entry (PS100001)

E. coli Selection: Kanamycin (25 ug/mL)

Fully Sequenced ORF: >MC211148 representing NM_001141983
Red=Cloning site Blue=ORF Orange=Stop codon

TTTGTGAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGGCCACTGAGGTCCACAATCTTCAGGAGCTCAGGCGAAGCGCTTCGCTGGCCACCAAGGTCTTTATCC
AGAGAGACTACAGTGACGGCACCCTGTGTCAGTTCCAGACCAAAATCCCGCCAGAGCTGGATAGCCGGAT
TGAGCGGCAGCTCTTTGAGGAGACGGTGAAGACCCTCAATGGTTTTATGCGGAGGCTGAGAAGATCGGG
GGCAGCTCCTATCTGGAGGGCTGCCTGGCCTGCGCCACGGCCTACTTCATCTTCCTCTGCATGGAGACCC
ACTATGAGAAGGTTCTCAAGAAGATTTCCGTTACATCCAGGAGCAGAATGAGAAGGTCTTTGCTCCCCG
AGGGCTCCTGCTCACGGATCCGGTGAACGCGGAATGAGGGTTATCGAGATCTCCATCTATGAGGACCGG
TGCAGCAGTGGTAGTTCCAGTAGCGGCAGCAGCAGCGGCAGCAGCAGCGCTGGCGGGGAGGGG
CGGGGGCCCGGTGA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI

ACCN: NM_001141983

Insert Size: 504 bp



OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	NM_001141983.1 , NP_001135455.1
RefSeq Size:	2839 bp
RefSeq ORF:	504 bp
Locus ID:	71146
UniProt ID:	Q9D428
Cytogenetics:	19 C3
Gene Summary:	<p>May be involved in protein transport from Golgi to cell surface.[UniProtKB/Swiss-Prot Function]</p> <p>Transcript Variant: This variant (1) represents the longer transcript. Both variants 1 and 2 encode the same protein.</p>