

Product datasheet for **MC211158**

Daw1 (NM_027725) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Daw1 (NM_027725) Mouse Untagged Clone
Tag: Tag Free
Symbol: Daw1
Synonyms: 4930563E19Rik; 4933429D11Rik; b2b1116Clo; b2b1584Clo; Wdr69
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >MC211158 representing NM_027725
Red=Cloning site **Blue**=ORF **Orange**=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGAAGCTTAAGAGCCTCTTGCTTCGCTATTACCCACCAGGAATTATGCTGGAATATGAAAAAGGTGGAG
AATTAAGACGAAATCCATAGATTTACTTGAGCTCAGTCCTAGTACTGATGTCAACACCCTAGTAGGAGA
GATCCAGCAAGCAGAGCCTCTGATCACAGCGTCCCGGACAAAGCAAGTCAGACTCCTAGTCAGAGGCTT
CAGGAGAACTCAGACAGCATAGCGACCACAATTCTACCTGTTAAGGTTCTCAGAGCACACATATTGC
CATTGACCAACGTGGCACTCAACAAAGCAGGCTCCTGCTTTATCACCGGAAGCTATGATCGGACCTGCAA
GGTGTGGGACACGGCGTCTGGAGAGGAGCTGCACACCCTGGAGGGCCACAAGAATGTGGTGTATGCCATA
GCCTTCAACAATCCTTACGGGAAGGTGCACACTTTAATTGGTCACTGTGCTGAGATCAGCAGCGCCTTGT
TCAATTGGGATTGCTCTAATACTAACAGGCTCCATGGACAAGACCTGCATGTTGTGGGATGCTACAAG
TGGGAAATATGTAGCAACCTTAACAGGCCATGACGATGAGATCCTAGACAGCTGCTTTGATTACTGGA
AACTTATCGCAACAGCTTCTGCCGATGGAACGGCAAGAGTTACAATGCAACCACAAGAAAGTGTGTA
CCAAACTGGAAGGCCATGAAGGAGAAATCTCAAAGATCTCCTTCAATCCCCAAGGAAACCGGCTCCTGAC
GGCAGCTCTGATAAGACAGCCAGAATCTGGGATGTTCAAGCTGGCCAGTGCCTCCAGGTCCTTGAGGGG
CACACTGATGAAATATTTTCATGTGCTTCAACTATAAAGGGAATATAGTCATTACAGGCAGCAAGGACA
ACTCCTGCAGGATATGGCGT**TGA**

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: Sgfl-MluI
ACCN: NM_027725
Insert Size: 933 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_027725.3 , NP_082001.2
RefSeq Size:	1371 bp
RefSeq ORF:	933 bp
Locus ID:	71227
Cytogenetics:	1 C5