

Product datasheet for MC211225

1700001F09Rik (NM_027940) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: 1700001F09Rik (NM_027940) Mouse Untagged Clone
Tag: Tag Free
Symbol: 1700001F09Rik
Synonyms: Gm10377
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >MC211225 representing NM_027940
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**GCGATCGCC**

ATGGGTTCCCAAGCAGGGTCATTTAGGAAGGCATCACCCAACTCCCAACATCGATGAGAACGAGAAAA
 GAATCAAGAGGTTGGAGAAGCTCAAAGAGATCTCCAGAACATTAATAAGAGAGACATTCTCCAGGG
 AATCCTGGCCAAGTACAAGGATTTGAATGACAGGACCAACTTTGAGACATTCATGCTTGAGATGCAACAT
 AACCAGATGATGACTGACCTGAAGAGAATGTCCCAGGACATCAGTGAGGCCTTGTACAAGTGCCAGCATT
 TGACTAAAGAAAATCAATTATACTGCTGCAGGAAGTCCACCTCCTGATTGAATCTAACCTCATCCAGCA
 CAAAGTCAGGATGTTGTGGAAGGAGAACAGACAGCTGCTAAGGAAGCAGATTGCATTGGAAGAATGCAAC
 ATAAAAACAAAAATACTATGTAAGAAGGCAGCCAGAAGATCAAAGACGAATATAGCAGGCAAGCAGCAG
 GAAAACAAATAATTGTGGTACCTGGAAGTAAAGAAGAGGCTGCTATGTCAAATCCATGTGCCT**G**A

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_027940
Insert Size: 555 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_027940.2, NP_082216.1</u>
RefSeq Size:	883 bp
RefSeq ORF:	555 bp
Locus ID:	71826
Cytogenetics:	14 B