

Product datasheet for MC211235

Gpr160 (NM_001134385) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Gpr160 (NM_001134385) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Gpr160
Synonyms:	1700025D19Rik; GPCR150
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>MC211235 representing NM_001134385 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGACAGCTCTTTCTTCAAAGAACTGCTCTTCAATACCAGCTCCATCAGTCTCCGAGCTCCTGGAGG
CCTCCTGTCTGTTCTTGATCATCCTGGGAAAGTGTATTAACATCCTCCTCTAAGAGTGAGAAG
GGGAGACGCCCGCTGGACTTTGATGGAGTACTTCTGCTTCTCACTCGCACTTGTGATCTCCTGCTGTTG
GTAACATTTCCATCCTGACCTACTTCAGGACTTTGTAGTTCTAGGTATTAGGTTTACCAGATACCACA
TCTGTCTGCTCACACAAATTATTTCTTCACTTATGGCTTCTTGCAATACCCAGTCTGCTCACTGGCTTG
TATAGACTACTGGTGAATCTGTCCAGAGCCAGCAAGCAGTCATCTAGGTGGCAGAAATTAATTTATTTT
TTGACAGTCATTTAACTTGATCTCAGTCCTTGCTTACGTTTTGGTTGACCCGGCCATCTCAGTCAGCC
TAAAGGCACACAGGGGCTATGTGTACCAGTGCCCGCCTATGTCAGCACTCAGAGTCACTGGCTCTCACT
GTCTATGCTGATGGTTTTATTTGGCTTTTCTGATTTTATGTCAGGAAGTTGTTGCCTTGCTTCAGGCT
ATGAGGATAGCATCCTATAAGAGCAAGGCCCTCTACTTTCTTTCCCACTCCACTGCGGTTACGCC
TGAGCTGCAGAGAAGCACTTTGCCAGGCTCATCGTGTGCTTTCTTGGACCTGGTTTCCCTTTGTAGC
ACTTCAGTCCCTCATCCTCTCCCTCAGAGTTTATGAGTCCCGGCTACATAGAAATGAACGTCCTGGCTG
TACTTTGTCAACAGTTTCTCATCGTGCAGTTTATTGGTTAACTGCCACAAGCTTGAATTAAGAGACA
GCTCGTTACCTGTGGATCCTTTCATCAACTGGAATGCTGCTTTGTTCCAGTCCATAGACTTAAGCAAGT
GGAGAGGCCCATGTCCATAGTCATCTGT**TAG**

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI



[View online »](#)

ACCN:	NM_001134385
Insert Size:	1011 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001134385.2</u> , <u>NP_001127857.1</u>
RefSeq Size:	2208 bp
RefSeq ORF:	1011 bp
Locus ID:	71862
UniProt ID:	<u>Q3U3F9</u>
Cytogenetics:	3 A3
Gene Summary:	Orphan receptor.[UniProtKB/Swiss-Prot Function]