

Product datasheet for MC211314

181004IL15Rik (NM_001163145) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids

Tag: Tag Free

Symbol: 181004IL15Rik

Synonyms: Kiaa1644; Shisall

Mammalian Cell Selection: Neomycin

Vector: pCMV6-Entry (PS100001)

E. coli Selection: Kanamycin (25 ug/mL)

Fully Sequenced ORF: >MC211314 representing NM_001163145
Red=Cloning site Blue=ORF Orange=Stop codon

TTTGTGAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGATGACCAGTTGTGGCCAGCGGTCCCGAACGTGCTCGCAGTATTCTCCCTGTTGTTTCTGCAGTCC
TGTCGCCACATTTCCGGGTCTGTGAGCCTTACACGGATCACAAAGGCCGCTACCACTTTGGCTTCCACTG
CCCCCGGCTCTCGGACAACAAACCTTCGTCTCTGCTGTCAACATAACAACACTGTCTTCAAATACTGC
TGCAACGAGACGGAGTTCCAGGCGGTATGCAGGCAAACTCACGGCCGGCCCCGAGGGCTACATGCACA
ACAACTACACGGCGCTTCTGGGAGTGTGATTTACGGCTTCTTCGTGCTCACGCTGCTGGTCTCTGGATCT
GCTCTATTATTCCGCCATGAACTACGATATTTGCAAGGTTTACCTCACTCGGTGGGGGATCCAGGGACGG
TGGATGAAACAGGACCCCGGCGGTGGGGGAACCCTGCTCGAGCCCTCGGCCAGGACAGCCGCCCCAC
AGCCGCAGCTCCCCCTGGTACCCTGCCGAAGCCCCCAGGCTGTGCACACTGCGGGGAGACACGCA
CAGCCACCCTGATGACCTTCCAGAGCTCATCCGCC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Chromatograms: https://cdn.origene.com/chromatograms/ja1803_e12.zip

Restriction Sites: SgfI-MluI

ACCN: NM_001163145

Insert Size: 600 bp



OTI Disclaimer:	Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery. The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	BC062953 , AAH62953
RefSeq Size:	5921 bp
RefSeq ORF:	600 bp
Locus ID:	72301
UniProt ID:	Q0VBP7
Cytogenetics:	15 E2